





## About DMS


Direct Medical Solution is a leading medical accessories manufacturer including SpO2, ECG/EKG, NIBP/IBP, Temp, EEG, ESU and Fetal. Established in 2008 with a registered capital of 10 million, headquartered in the coastal city Shenzhen, which is regarded as the 'Silicon Valley' of the China's medical electronics industry.

DMS has a complete system of R&D, Production and Sales with over decades years development. We have been adhering to the 'innovation & responsibility' as our company philosophy and dedicate our selves to an one-stop solution for medical accessories with reliable quality and wide compatibility to our customers. Meanwhile, DMS approved various certifications like TUV, CE&FDA, which allow us access to the major global markets, and helps our distributors to succeed in more than 100 countries and regions around the world.

## | Directory

<a href="#">SpO2 Sensors &amp; Adapter Cables</a>	001
<a href="#">ECG/EKG Cables &amp; Accessories</a>	039
<a href="#">NIBP Cuffs &amp; Accessories</a>	080
<a href="#">IBP Cable &amp; Transducer</a>	097
<a href="#">Temperature Probes</a>	104
<a href="#">Fetal Probes</a>	115
<a href="#">ESU Accessories</a>	119
<a href="#">Veterinary Accessories</a>	123
<a href="#">Respiratory Consumables</a>	125





## Disposable SpO2 Sensors

- Adhesive Tape and Foam are available
- Latex-free
- Available in adult, pediatric, infant and Neonate sizes
- Single patient use to reduce cross-contamination



## Disposable SpO2 Sensors


### Certification

Free Sales Certificate(FSC) No.:SAMD2021-046

### Compliance with

- ISO 10993-1:2018
- ISO 10993-5:2009
- ISO 10993-10:2010
- IEC 62366-1:2015
- EN 60601-1:2006/A1:2013
- EN ISO 80601-2-56:2017/A1:2020

Product types	Photo	P/N Example	P/N Guide
Adhesive Tape sensor		DMS 1019	P: SpO2 1: disposable sensors 0: Neonate/adult adhesive sensor (1: Infant; 2: Pediatric; 3: Adult) 19: Nellcor non-oximax 19A: Nellcor oximax
Non-Adhesive Foam sensor		DMS 1819	P: SpO2 1: disposable sensors 8: Neonate/adult foam non-adhesive sensor (9: Infant/Pediatric ) 19: Nellcor non-oximax 19A: Nellcor oximax
Adhesive Foam sensor		DMS 1419	P: SpO2 1: disposable sensors 4: Neonate/adult foam adhesive sensor (5: Infant; 6: Pediatric; 7: Adult ) 19: Nellcor non-oximax 19A: Nellcor oximax
Warranty: no warranty for disposable sensors except defectiveness found in the first use Package: 1pcs/bag, 25pcs/box			



### Disposable SpO2 Sensors

- Adhesive Tape and Foam are available
- Latex-free
- Available in adult, pediatric, infant and Neonate sizes
- Single patient use to reduce cross-contamination



## Disposable SpO2 Sensors

Product types	Photo	P/N Example	P/N Guide
Adhesive Tape sensor		DMS 1019	P: SpO2 1: disposable sensors 0: Neonate/adult adhesive sensor (1: Infant; 2: Pediatric; 3: Adult) 19: Nellcor non-oximax 19A: Nellcor oximax
Non-Adhesive Foam sensor		DMS 1819	P: SpO2 1: disposable sensors 8: Neonate/adult foam non-adhesive sensor (9: Infant/Pediatric ) 19: Nellcor non-oximax 19A: Nellcor oximax
Adhesive Foam sensor		DMS 1419	P: SpO2 1: disposable sensors 4: Neonate/adult foam adhesive sensor (5: Infant; 6: Pediatric; 7: Adult ) 19: Nellcor non-oximax 19A: Nellcor oximax

Warranty: no warranty for disposable sensors except defectiveness found in the first use  
 Package: 1pcs/bag, 25pcs/box

Disposable SpO2 Sensors

Monitor	Adhesive Tape	Adhesive Foam	Non-Adhesive Foam	Size	Len.	Conn. P/N
BCI	DMS1003 1302	DMS1403	DMS1803	Neonate/Adult 成人/新生儿 ( < 3kg/>40kg)	0.9m	 CN81
	DMS1103 1303	DMS1503	DMS1903	Infant 婴儿 (3-20kg)		
	DMS1203 1301	DMS1603	DMS1903	Pediatric 儿童 (10-50kg)		
	DMS1303 1300	DMS1703	DMS1803	Adult 成人 (>30kg)		
CSI	DMS1007 573SD	DMS1407	DMS1807	Neonate/Adult 成人/新生儿 ( < 3kg/>40kg)	0.9m	 CN06
	DMS1107 572SD	DMS1507	DMS1907	Infant 婴儿 (3-20kg)		
	DMS1207 571SD	DMS1607	DMS1907	Pediatric 儿童 (10-50kg)		
	DMS1307 570SD	DMS1707	DMS1807	Adult 成人 (>30kg)		
Datascope	DMS1008 0998-00-0074-04	DMS1408	DMS1808	Neonate/Adult 成人/新生儿 ( < 3kg/>40kg)	0.9m	 CN08
	DMS1108 0998-00-0074-03	DMS1508	DMS1908	Infant 婴儿 (3-20kg)		
	DMS1208 0998-00-0074-04	DMS1608	DMS1908	Pediatric 儿童 (10-50kg)		
	DMS1308 0998-00-0076-05	DMS1708	DMS1808	Adult 成人 (>30kg)		
GE-Datex	DMS1010 SAS-AF	DMS1410	DMS1810	Neonate/Adult 成人/新生儿 ( < 3kg/>40kg)	0.9m	 CN04
	DMS1110 SAS-AF	DMS1510	DMS1910	Infant 婴儿 (3-20kg)		
	DMS1210 SAS-AF/SAS-AF	DMS1610	DMS1910	Pediatric 儿童 (10-50kg)		
	DMS1310 SAS-AF/SAS-AF	DMS1710	DMS1810	Adult 成人 (>30kg)		
GE-OXYTIP+	DMS1010L OXY-AF-10	DMS1410L	DMS1810L	Neonate/Adult 成人/新生儿 ( < 3kg/>40kg)	0.9m	 CN17
	DMS1110L OXY-AF-10	DMS1510L	DMS1910L	Infant 婴儿 (3-20kg)		
	DMS1210L OXY-AF-10	DMS1610L	DMS1910L	Pediatric 儿童 (10-50kg)		
	DMS1310L OXY-AF-10	DMS1710L	DMS1810L	Adult 成人 (>30kg)		

Disposable SpO2 Sensors

Monitor	Adhesive Tape	Adhesive Foam	Non-Adhesive Foam	Size	Len.	Conn. P/N
GE-TruSignal	DMS1010S TS-AF-10 TS-AF-25	DMS1410S	DMS1810S	Neonate/Adult ( < 3kg/>40kg)	0.9m	 CN64
	DMS1110S TS-AF-10 TS-AF-25	DMS1510S	DMS1910S	Infant (3-20kg)		
	DMS1210S TS-AF-10 TS-AF-25	DMS1610S	DMS1910S	Pediatric (10-50kg)		
	DMS1310S TS-AF-10 TS-AF-25	DMS1710S	DMS1810S	Adult (>30kg)		
Masim	DMS1015A 1862	DMS1415A	DMS1815A	Neonate/Adult ( < 3kg/>40kg)	0.9m	 CN23
	DMS1115A 1861	DMS1515A	DMS1915A	Infant (3-20kg)		
	DMS1215A 1860	DMS1615A	DMS1915A	Pediatric (10-50kg)		
	DMS1315A 1859	DMS1715A	DMS1815A	Adult (>30kg)		
	DMS1015A	DMS1415A	DMS1815A	Neonate/Adult ( < 3kg/>40kg)		 CN23
	DMS1115A 2329	DMS1515F	DMS1915F	Infant (3-20kg)		
	DMS1215A	DMS1615F	DMS1915F	Pediatric (10-50kg)		
	DMS1315A	DMS1715F	DMS1815F	Adult (>30kg)		
	DMS1015B 1002	DMS1415B	DMS1815B	Neonate/Adult ( < 3kg/>40kg)		 CN24
	DMS1115B	DMS1515B	DMS1915B	Infant (3-20kg)		
	DMS1215B 1025	DMS1615B	DMS1915B	Pediatric (10-50kg)		
	DMS1315B 1001	DMS1715B	DMS1815B	Adult (>30kg)		
	DMS1015S 4003	DMS1415S	DMS1815S	Neonate/Adult ( < 3kg/>40kg)		 CN443
	DMS1115S	DMS1515S	DMS1915S	Infant (3-20kg)		
	DMS1215S	DMS1615S	DMS1915S	Pediatric (10-50kg)		
	DMS1315S	DMS1815S	DMS1815S	Adult (>30kg)		

Disposable SpO2 Sensors

Monitor	Adhesive Tape	Adhesive Foam	Non-Adhesive Foam	Size	Len.	Conn. P/N
Masim	DMS1015H	DMS1415H	DMS1815H	Neonate/Adult (< 3kg/>40kg)	0.9m	 CN86
	DMS1115H	DMS1515H	DMS1915H	Infant (3-20kg)		
	DMS1215H	DMS1615H	DMS1915H	Pediatric (10-50kg)		
	DMS1315H	DMS1715H	DMS1815H	Adult (>30kg)		
Nellcor	DMS1019 N25	DMS1419	DMS1819	Neonate/Adult (< 3kg/>40kg)	0.9m	 CN30
	DMS1119 I20	DMS1519	DMS1919	Infant (3-20kg)		
	DMS1219 D20	DMS1619	DMS1919	Pediatric (10-50kg)		
	DMS1319 D25	DMS1719	DMS1819	Adult (>30kg)		
	DMS1019A MAX-N	DMS1419A	DMS1819A	Neonate/Adult (< 3kg/>40kg)	0.9m	 CN31
	DMS1119A MAX-I	DMS1517A	DMS1919A	Infant (3-20kg)		
	DMS1219A MAX-P	DMS1619A	DMS1919A	Pediatric (10-50kg)		
	DMS1319A MAX-A	DMS1719A	DMS1819A	Adult (>30kg)		
Nihon Kohden	DMS1021 TL-253T	DMS1421	DMS1821	Neonate/Adult (< 3kg/>40kg)	0.9m	 CN33
	DMS1121 TL-252T	DMS1521	DMS1921	Infant (3-20kg)		
	DMS1221 TL-252T	DMS1621	DMS1921	Pediatric (10-50kg)		
	DMS1321 TL-251T	DMS1721	DMS1821	Adult (>30kg)		
Nonin	DMS1022 9000N	DMS1422	DMS1822	Neonate/Adult (< 3kg/>40kg)	0.9m	 CN83
	DMS1122 9000I	DMS1522	DMS1922	Infant (3-20kg)		
	DMS1222 9000P	DMS1622	DMS1922	Pediatric (10-50kg)		
	DMS1322 9000A	DMS1722	DMS1822	Adult (>30kg)		
Novamatrix	DMS1023 AS140	DMS1423	DMS1823	Neonate/Adult (< 3kg/>40kg)	0.9m	 CN30
	DMS1123 AS120	DMS1523	DMS1923	Infant (3-20kg)		
	DMS1223 AS120	DMS1623	DMS1923	Pediatric (10-50kg)		
	DMS1323 AS110	DMS1723	DMS1823	Adult (>30kg)		



Reusable SpO2 Sensors

# Reusable SpO2 Sensors

- Finger clip, Soft tip, Wrap, Multi-Y, Ear clip types available
- Latex-free
- Available in adult, pediatric, infant and neonate sizes
- Different types are Compatible with most of patient monitors
- High performance/price ratio for Reusable SpO2 sensor to save cost
- Easy to clean and repair



Product types	Photo	P/N Example	P/N Guide
Reusable SpO2 Sensors		DMS8325A	P: SpO2 8: Adult soft (9: Adult clip; 7: Pediatric clip; 6: Pediatric soft; 5: Neonate wrap; 4: Multi Y; 3: Ear clip) 3: 3.0m (1: 1.0m) 25A: Philips 8 pins (25: Philips 12 pins; 19: Nellcor 7 pins; 19A: Nellcor 9 pins)
Warranty: one-year warranty for reusable sensors and cables Package: 1pcs/bag			

# Reusable SpO2 Sensors

Monitor	Adult clip	Adult soft	Pedia. Clip	Pedia. Soft	Neonate wrap	Multi-Y	Adult Ear clip	Connector/Adapter/Len.
<b>Artemstrong Medica</b>								
AD-900, 1000, 2000, 3000, 4000(BCI tech)	DMS9303A	DMS8303A	DMS7303A	DMS6303A	DMS5303A	DMS4303A	DMS3303A	CN81, Direct ,3m
<b>Artema</b>								
Diascope	DMS9301A	DMS8301A	DMS7301A	DMS6301A	DMS5301A	DMS4301A	DMS3301A	CN01, Direct ,3m
S&W	DMS9301A SI033057	DMS8301A	DMS7301A	DMS6301A	DMS5301A	DMS4301A	DMS3301A	CN02, Direct ,3m
<b>Baxter/Simed</b>								
ASAT	DMS9302A 600015	DMS8302A	DMS7302A	DMS6302A	DMS5302A	DMS4302A	DMS3302A	CN03, Direct ,3m
S-50, S100, S-100e	DMS9302	DMS8302	DMS7302	DMS6302	DMS5302	DMS4302	DMS3302	CN03, Direct ,3m
Monarch/Explorer	DMS9319	DMS8319	DMS7319	DMS6319	DMS5319	DMS4319	DMS3319	CN30, Direct ,3m
<b>BCI/Smiths</b>								
3100, 3101, 6100, 9100, Advisor, Autocorr, Mini Torr, Mini Torr Plus	DMS9303	DMS8303	DMS7303	DMS6303	DMS5303	DMS4303	DMS3303	CP7, Direct ,3m
3300, 3301, 3302, 3303, 3304, 3401, 3180 ChoiceMed MD2000A	DMS9103A	DMS8103A	DMS7103A	DMS6103A	DMS5103A	DMS4103A	DMS3103A	CN81, Direct ,1m
	DMS9303A 3044/3444	DMS8303A	DMS7303A	DMS6303A	DMS5303A 3025/3026	DMS4303A	DMS3303A 3043	CN81, Direct ,3m
<b>Biosys</b>								
BPM-200, Vital Sign Monitor SENTRY	DMS9304	DMS8304	DMS7304	DMS6304	DMS5304	DMS4304	DMS3304	CN05, Direct ,3m
<b>Bird</b>								
4400	DMS9307	DMS8307	DMS7307	DMS6307	DMS5307	DMS4307	DMS3307	CM5, Direct ,3m
<b>Biocare</b>								
Biocare iM12 (BCI)	DMS9393	DMS8393	DMS7393	DMS6393	DMS5393	DMS4393	DMS3393	CN95, Direct ,3m



Reusable SpO2 Sensors

Monitor	Adult clip	Adult soft	Pedia. Clip	Pedia. Soft	Neonate wrap	Multi-Y	Adult Ear clip	Connector/Adapter/Len.
<b>Bionet</b>								
BM3(before 2005), BM3 Vet	DMS9306	DMS8306	DMS7306	DMS6306	DMS5306	DMS4306	DMS3306	 CP7-1, Direct ,3m
Bm3	DMS9306A	DMS8306A	DMS7306A	DMS6306A	DMS5306A	DMS4306A	DMS3306A	 CN63, Direct ,3m
BM5	DMS9306B	DMS8306B	DMS7306B	DMS6306B	DMS5306B	DMS4306B	DMS3306B	 CN78, Direct ,3m
<b>Biolight</b>								
Biolight M12; LGMD (BCI)	DMS9305	DMS8305	DMS7305	DMS6305	DMS5306	DMS4305	DMS3305	 CPS, Direct ,3m
Biolight M69; M9000 Biocare Neusoft NSC-M10; Contec CMS9000, CMS6000 (Digital)	DMS9305A	DMS8305A	DMS7305A	DMS6305A	DMS5306A	DMS4305A	DMS3305A	
M700,M800, Creative SP20 (Digital)	DMS9105C	DMS8105C	DMS7105C	DMS6105C	DMS5105C	DMS4105C	DMS3105C	 CN30, P0205, 1m
M7000 ,smk-1 , Advanced AM1500 (Digital)	DMS9305D	DMS8305D	DMS7305D	DMS6305D	DMS5305D	DMS4305D	DMS3305D	 CN61, Direct ,3m
Schiller,Argus,LSM (Nellcor Oximax)	DMS9305F	DMS8305F	DMS7305F	DMS6305F	DMS5305F	DMS4305F	DMS3305F	
Analogic medical AM 1500 (Digital new) (N)	DMS9305N	DMS8305N	DMS7305N	DMS6305N	DMS5305N	DMS4305N	DMS3305N	
(Digital)	DMS9305E	DMS8305E	DMS7305E	DMS6305E	DMS5305E	DMS4305E	DMS3305E	 CN62, Direct ,3m
A8,M9000,V6,Advanced PM-200M,Any View A5 (Digital)	DMS9305H	DMS8305H	DMS7305H	DMS6305H	DMS5305H	DMS4305H	DMS3305H	 CN82, Direct ,3m
A8 ,M9000,V6, Advanced PM-200M, Any View A5,Innomed T12Plus (Digital new) (N)	DMS9305J	DMS8305J	DMS7305J	DMS6305J	DMS5305J	DMS4305J	DMS3305J	
A8 (Nellcor Oximax)	DMS9305G	DMS8305G	DMS7305G	DMS6305G	DMS5305G	DMS4305G	DMS3305G	 CN30-C, Direct ,1m
AM1500 (Digital new) (N)	DMS9105K	DMS8105K	DMS7105K	DMS6105K	DMS5105K	DMS4105K	DMS3105K	
Anyview (Analog new) (N)	DMS9105M	DMS8105M	DMS7105M	DMS6105M	DMS5105M	DMS4105M	DMS3105M	

Reusable SpO2 Sensors

Monitor	Adult clip	Adult soft	Pedia. Clip	Pedia. Soft	Neonate wrap	Multi-Y	Adult Ear clip	Connector/Adapter/Len.
<b>Biolight</b>								
Anyview S12 (Analog new) (N)	DMS9305L	DMS8305L	DMS7305L	DMS6305L	DMS5305L	DMS4305L	DMS3305L	 CN82B, Direct ,3m
<b>CAS Medical</b>								
OscilloMate 9002, 9303, 9302S (Nonin tech)	DMS922B	DMS822B	DMS722B	DMS6322B	DMS5322B	DMS4322B	DMS3322B	 CN83, Direct ,3m
740-750 (Masim tech)	DMS9315	DMS8315	DMS7315	DMS6315	DMS5315	DMS4315	DMS3315	 CN22, Direct ,3m
<b>CSI</b>								
501, 503, 504, 507, 508, 581, 1100,8100, Poet(602-3A), Scholar	DMS9307 934-10LN	DMS8307 975AD-10L	DMS7307 1563-10L	DMS6307 975PD-10L	DMS5307	DMS4307 940-10L	DMS3307 1567-10L	 CM5, Direct ,3m
90 degree angle	DMS9107A 934SDN	DMS8107A 975AD-3	DMS7107A 1563-SD	DMS6107A 975PD-3	DMS5107A	DMS4107A 940SD	DMS3107A 1567SD	 CN06, See Cable, 1m
503 Spot, 504Plus, 504DX, 506DXN 90 degree angle	DMS9307A 934-10DN	DMS8307A 975AD-10	DMS7307A 1563-10D	DMS6307A 975PD-10	DMS5307A	DMS4307A 940-10D	DMS3307A 1567-10D	 CN06, Direct ,3m
<b>Colin</b>								
BP88/306/508/608/BX10ne (Nellcor tech)	DMS9319	DMS8319	DMS7319	DMS6319	DMS5319	DMS4319	DMS3319	 CN30, Direct ,3m
BP 88S	DMS9315	DMS8315	DMS7315	DMS6315	DMS5315	DMS4315	DMS3315	 CN22, Direct ,3m
Bx10ma (Masim tech)	DMS9319B	DMS8319B	DMS7319B	DMS6319B	DMS5319B	DMS4319B	DMS3319B	
<b>Comen</b>								
Star8000C/8000D , Star 8000 ,C50 (Digital)	P9361A	P8361A	P7361A	P6361A	P5361A	P4361A	P3361A	 CP6, Direct ,3m
CM8000A,Star-1000(BCI)	DMS9361	DMS8361	DMS7361	DMS6361	DMS5361	DMS4361	DMS3361	
DOTOP, EDANINS EM9000E, SPRINGS SPR9000G, Jumbor, Huateng H8, AVEUS AV-Pro(BCI)	DMS9361B	DMS8361B	DMS7361B	DMS6361B	DMS5361B	DMS4361B	DMS3361B	 CP6R, Direct ,3m
(Digital)	DMS9361C	DMS8361C	DMS7361C	DMS6361C	DMS5361C	DMS4361C	DMS3361C	

Reusable SpO2 Sensors

Monitor	Adult clip	Adult soft	Pedia. Clip	Pedia. Soft	Neonate wrap	Multi-Y	Adult Ear clip	Connector/Adapter/Len.
<b>Comen</b>								
C60 (Nellcor Oximax)	DMS9361D	DMS8361D	DMS7361D	DMS6361D	DMS5361D	DMS4361D	DMS3361D	 CP8R, Direct ,3m
C60 (Masim tech)	DMS9361E	DMS8361E	DMS7361E	DMS6361E	DMS5361E	DMS4361E	DMS3361E	
STAR 8000A (BCI tech)	DMS9361J	DMS8361J	DMS7361J	DMS6361J	DMS5361J	DMS4361J	DMS3361J	 CP6S, Direct ,3m
C60 (Masim tech)	DMS9361K	DMS8361K	DMS7361K	DMS6361K	DMS5361K	DMS4361K	DMS3361K	 CN22, Direct ,3m
C60 (Masim tech)	DMS9361L	DMS8361L	DMS7361L	DMS6361L	DMS5361L	DMS4361L	DMS3361L	 CM8-1B, Direct ,3m
C50/C60/C80(Masim tech with M* chips)	DMS9361G	DMS8361G	DMS7361G	DMS6361G	DMS5361G	DMS4361G	DMS3361G	 CN93, Direct ,3m
Masim	DMS9361P	DMS8361P	DMS7361P	DMS6361P	DMS5361P	DMS4361P	DMS3361P	
C100A Telemetry Monitoring (Digital)	DMS9365D	DMS8365D	DMS7365D	DMS6365D	DMS5365D	DMS4365D	DMS3365D	 CN96, Direct ,3m
C30/C80(Nellcor tech), Star8000 NCS(Nellcor tech)	DMS9361H	DMS8361H	DMS7361H	DMS6361H	DMS5361H	DMS4361H	DMS3361H	
<b>Contec</b>								
CMS7000/ 8000/9000 Sinohero H8	DMS93132	DMS83132	DMS73132	DMS63132	DMS53132	DMS43132	DMS33132	 CN452, Direct ,3m
<b>Choice Electronic</b>								
All (BCI tech)	DMS9303A	DMS8303A	DMS7303A	DMS6303A	DMS5303A	DMS4303A	DMS3303A	 CN81, Direct ,3m
M3,MPC300, MMED6000DP	DMS9342	DMS8342	DMS7342	DMS6342	DMS5342	DMS4342	DMS3342	 CP6, Direct ,3m
Ecomed EPM-50, EPM-80 iplus (Digital)	DMS9342A	DMS8342A	DMS7342A	DMS6342A	DMS5342A	DMS4342A	DMS3342A	
MD2000B	DMS9342C	DMS8342C	DMS7342C	DMS6342C	DMS5342C	DMS4342C	DMS3342C	 CN99, Direct ,3m








Reusable SpO2 Sensors

Monitor	Adult clip	Adult soft	Pedia. Clip	Pedia. Soft	Neonate wrap	Multi-Y	Adult Ear clip	Connector/Adapter/Len.
<b>Charmcare</b>								
Charmcare Cx100 (Nellcor tech)	DMS9319	DMS8319	DMS7319	DMS6319	DMS5319	DMS4319	DMS319	 CN30, Direct ,3m
<b>Cardioline</b>								
Cardioline Vsign 200B(BCI)	DMS9354	DMS8354	DMS7354	DMS6354	DMS5354	DMS4354	DMS3354	 CP6, Direct ,3m
<b>Creative</b>								
Creative, Heal force: advance 120, Deluxe 70	DMS9395	DMS8395	DMS7395	DMS6395	DMS5395	DMS4395	DMS3395	 CP8, Direct ,3m
<b>CU Medical</b>								
LiFEGAIN CU-HD1	DMS9103	DMS83103	DMS73103	DMS63103	DMS53103	DMS43103	DMS33103	 CP8, Direct ,3m
<b>Dolphin</b>								
2100, 2150 (Digital)	DMS9162A	DMS8162A	DMS7162A	DMS6162A	DMS5162A	DMS4162A	DMS3162A	 CN71, See Cable, 1m
	DMS9362	DMS8362	DMS7362	DMS6362	DMS5362	DMS4362	DMS3362	 CN22, Direct ,3m
<b>Drager-Siemens</b>								
Drager Oxipac 2000, 2500 (BCI tech)	DMS9303A	DMS8303A	DMS7303A	DMS6303A	DMS5303A	DMS4303A	DMS3303A	 CN81, Direct ,3m
Vitalert 1000, Cicero, GammaXL, DeltaXL, Kappa	DMS9309	DMS8309	DMS7309	DMS6309	DMS5309	DMS4309	DMS3309	 CM8, Direct ,3m
Sirecust, Sirem Spo2	DMS9309A	DMS8309A	DMS7309A	DMS6309A	DMS5309A	DMS4309A	DMS3309A	 CN09, Direct ,3m
Sc6000, 7000,9000	DMS9309B	DMS8309B	DMS7309B	DMS6309B	DMS5309B	DMS4309B	DMS3309B	 CN10, Direct ,3m
<b>DST(USA)</b>								
DST(USA)	DMS9381	DMS8381	DMS7381	DMS6381	DMS5381	DMS4381	DMS3381	 CP6R, Direct ,3m
<b>GE-Corometrics</b>								
Corometrics 188,129, 129F,128, 128F	DMS9319	DMS8319	DMS7319	DMS6319	DMS5319	DMS4319	DMS3319	 CN30, Direct ,3m
Corometrics 120, 510, 511, 118,556	DMS9310	DMS8310	DMS7310	DMS6310	DMS5310	DMS4310	DMS3310	 CN11, Direct ,3m












Reusable SpO2 Sensors

Monitor	Adult clip	Adult soft	Pedia. Clip	Pedia. Soft	Neonate wrap	Multi-Y	Adult Ear clip	Connector/Adapter/Len.
<b>GE-Critikon</b>								
Dinamap 8700/8710/8720, Plus, Oxyshuttle	DMS9310A	DMS8310A	DMS7310A	DMS6310A	DMS5310A	DMS4310A	DMS3310A	 CM7, Direct, 3m
Dinamap 9700 / 9710 / 9720, Plus, Compact, MPS, Select	DMS9319	DMS8319	DMS7319	DMS6319	DMS5319	DMS4319	DMS3319	 CN30, Direct, 3m
Dinamap ProSeries V1 and Pro 1000 V2 (Nellcor tech)	DMS9319B	DMS8319B	DMS7319B	DMS6319B	DMS5319B	DMS4319B	DMS3319B	 CN22, Direct, 3m
ProCare; ProSeries V2; Pro 1000 V3 (Masim)	DMS9315	DMS8315	DMS7315	DMS6315	DMS5315	DMS4315	DMS3315	
Same as above (Nellcor Oximax)	DMS9119A	DMS8119A	DMS7119A	DMS6119A	DMS5119A	DMS4119A	DMS3119A	 CN31, P0319C, 1m
<b>GE-Datex Ohmeda</b>								
S/5, AS/3, CS/3 Modular Monitors; with M-series modules; Cardiocap/5, Cardiocap S/5 Light Monitor; Capnomac Ultima Satlite Plus II	DMS9310B OXY-F4-N	DMS8310B	DMS7310B	DMS6310B	DMS5310B	DMS4310B	DMS3310B OXY-E4-N	 CN12, Direct, 3m
	DMS9110C OXY-F-DB	DMS8110C	DMS7110C	DMS6110C	DMS5110C OXY-W-DB	DMS4110C	DMS3110C OXY-E-DB	
M-NSAT Modules for AS/3, CS/3, S/5 Modular	DMS9319B	DMS8319B	DMS7319B	DMS6319B	DMS5319B	DMS4319B	DMS3319B	 CN22, Direct, 3m
E-NSAT Modules for i4 Modular Monitors	DMS9119A	DMS8119A	DMS7119A	DMS6119A	DMS5119A	DMS4119A	DMS3119A	 CN31, P0319C, 1m
<b>GE-Ohmeda 欧美达</b>								
Ohmeda TuffSat; 3775, 37xx Series, 3800, 3900, 3900P, 4700 and 5250 respiratory gas monitors; Cardiocap/5	DMS9310H OXY-F4-H	DMS8310H	DMS7310H	DMS6310H	DMS5310H OXY-W4-H	DMS4310H	DMS3310H OXY-E4-H	 CN16, Direct, 3m
Ohmeda TuffSat	DMS9110H OXY-F1-H	DMS8110H	DMS7110H	DMS6110H	DMS5110H	DMS4110H	DMS3110H	 CN16, Direct, 1m
Ohmeda TuffSat; 3775, 37xx Series, 3800, 3900, 3900P, 4700 and 5250 respiratory gas monitors; Cardiocap/5	DMS9110J	DMS8110J	DMS7110J	DMS6110J	DMS5110J	DMS4110J	DMS3110J	 CN04, P0210H, 1m
Ohmeda TruSat, Ohmeda TruSat TD	DMS9310M OXY-F4-MC	DMS8310M	DMS7310M	DMS6310M	DMS5310M	DMS4310M	DMS3310M OXY-E4-MC	 CN18, Direct, 3m
GE-Ohmeda OxyTip sensor	DMS9110L OXY-F-UN	DMS8110L	DMS7110L	DMS6110L	DMS5110L OXY-W-UN	DMS4110L OXY-SE-3	DMS3110L	 CN17, GE Cable, 1m




Reusable SpO2 Sensors

Monitor	Adult clip	Adult soft	Pedia. Clip	Pedia. Soft	Neonate wrap	Multi-Y	Adult Ear clip	Connector/Adapter/Len.
<b>GE-Ohmeda</b>								
S/5, i4 Modular Monitors with E-series modules; Corometrics 259 with TruSignal(Ohmeda tech) B20/30/40 Series	DMS9310K TS-F4-GE	DMS8310K	DMS7310K	DMS6310K	DMS5310K	DMS4310K	DMS3310K TS-E4-GE	 CN15, Direct, 3m
<b>GE-Marquette</b>								
Dash 2000/3000/4000; Solar8000/8000M,i/9500 with Tram Module451/451N/851/851N;Cardiolab/Mac-Labwith Tram Module 451/451N;Corometrics120/2120 series; Eagle 1000/3000/4000;Responder 3000 (Nellcor tech) Oxismart	DMS9310E	DMS8310E	DMS7310E	DMS6310E	DMS5310E	DMS4310E	DMS3310E	 CN14, Direct, 3m
Dash 3000/4000/5000 (OxiMax);Solar 8000/8000M, i9500 with Tram module 451N5/851N5	DMS9310F	DMS8310F	DMS7310F	DMS6310F	DMS5310F	DMS4310F	DMS3310F	 CN15, Direct, 3m
Dash 3000/4000/5000 Solar 8000/8000M, i/9500 with Tram Module 451M/851M Cardiolab /Mac-Lab with Tram Module 451M; Corometrics 259(Masim)	DMS9310G	DMS8310G	DMS7310G	DMS6310G	DMS5310G	DMS4310G	DMS3310G	
GE CARESCAPE V100, GE CARESCAPE B650 Masim Tech (with M* chip)	DMS9310W	DMS8310W	DMS7310W	DMS6310W	DMS5310W	DMS4310W	DMS3310W	
<b>GE-Hellige</b>								
Eagle 1000, 3000, SMU, SMK, SMV, cardioserv	DMS9310D	DMS8310D	DMS7310D	DMS6310D	DMS5310D	DMS4310D	DMS3310D	 CN13, Direct, 3m
<b>GE TruSignal</b>								
TruSignal sensor	DMS9110S TS-F-D	DMS8110S TS-SA-D	DMS7110S	DMS6110S TS-SP-D	DMS5110S TS-W-D	DMS4110S TS-SE-3	DMS3110S	 CN64, GE Cable, 1m
	DMS9110T	DMS8110T	DMS7110T	DMS6110T	DMS5110T	DMS4110T	DMS3110T TS-E-D	 CN74, GE Cable, 1m
<b>Guoteng</b>								
Carewell, Nova-3D (BCI tech)	DMS9350	DMS8350	DMS7350	DMS6350	DMS5350	DMS4110T	DMS3110T	 CP6, Direct, 3m












Reusable SpO2 Sensors

Monitor	Adult clip	Adult soft	Pedia. Clip	Pedia. Soft	Neonate wrap	Multi-Y	Adult Ear clip	Connector/Adapter/Len.
<b>Goldway</b>								
940 (Nellcor tech)	DMS9311	DMS8311	DMS7311	DMS6311	DMS5311	DMS4311	DMS3311	 CMS, Direct ,3m
Nanjing SD, SENMIND 900F, KELIWEI, STRONG 7000C, HRD EM6/EM7, Snolife Smart wave15, CMS6500	DMS9311A	DMS8311A	DMS7311A	DMS6311A	DMS5311A	DMS4311A	DMS3311A	 CP5, Direct ,3m
UT4000F (Digital)	DMS9311D	DMS8311D	DMS7311D	DMS6311D	DMS5311D	DMS4311D	DMS3311D	
(Nellcor tech)	DMS9311B	DMS8311B	DMS7311B	DMS6311B	DMS5311B	DMS4311B	DMS3311B	 CP7, Direct ,3m
UT4000B (Nellcor oximax)	DMS9311E	DMS8311E	DMS7311E	DMS6311E	DMS5311E	DMS4311E	DMS3311E	
<b>Hwatime</b>								
Hwatime UN-Medical galaxy	DMS9364A	DMS8364A	DMS7364A	DMS6364A	DMS5364A	DMS4364A	DMS3364A	 CP6R, Direct ,3m
H6(BCI)	DMS93126	DMS83126	DMS73126	DMS63126	DMS53126	DMS43126	DMS33126	
<b>Ruibo (Hunan)</b>								
16#12.5mm (BCI tech)	DMS9367	DMS8367	DMS7367	DMS6367	DMS5367	DMS4367	DMS3367	 CN68, Direct ,3m
12#9mm, PM-9000A+ (BCI tech)	DMS9367A	DMS8367A	DMS7367A	DMS6367A	DMS5367A	DMS4367A	DMS3367A	 CN73, Direct ,3m
12#9mm, PM-9000A+ (Digital)	DMS9367B	DMS8367B	DMS7367B	DMS6367B	DMS5367B	DMS4367B	DMS3367B	 CN420, Direct ,3m
<b>Infinium</b>								
OMNI 2, OMNI 3; STRONG(Digital)	DMS9365	DMS8365	DMS7365	DMS6365	DMS5365	DMS4365	DMS3365	 CP5R, Direct ,3m
OMNI(BCI Tech)	DMS9365A	DMS8365A	DMS7365A	DMS6365A	DMS5365A	DMS4365A	DMS3365A	 CP6, Direct ,3m
OMNI III(Digital)	DMS9365E	DMS8365E	DMS7365E	DMS6365E	DMS5365E	DMS4365E	DMS3365E	 CP5, Direct ,3m
OMNI II	DMS9365B	DMS8365B	DMS7365B	DMS6365B	DMS5365B	DMS4365B	DMS3365B	 CN90, Direct ,3m









Reusable SpO2 Sensors

Monitor	Adult clip	Adult soft	Pedia. Clip	Pedia. Soft	Neonate wrap	Multi-Y	Adult Ear clip	Connector/Adapter/Len.
<b>Infinium</b>								
OMNI III , Comen C50/80 (Digital)	DMS9365C	DMS8365C	DMS7365C	DMS6365C	DMS5365C	DMS4365C	DMS3365C	 CN93, Direct ,3m
<b>Ivy Biomedical</b>								
405P, 405T	DMS9319	DMS8319	DMS7319	DMS6319	DMS5319	DMS4319	DMS3319	 CN30, Direct ,3m
Vital-Guard 450C	DMS9315	DMS8315	DMS7315	DMS6315	DMS5315	DMS4315	DMS3315	 CN22, Direct ,3m
<b>Innomed</b>								
Oxycard, Innocare-T0	DMS9369	DMS8369	DMS7369	DMS6369	DMS5369	DMS4369	DMS3369	 CP7R, Direct ,3m
<b>Invivo</b>								
4500, 4500 Plus, Millenia Omega 1445,Omni-Trak TVS, Sat Scout, 4500 Plus 2	DMS9312	DMS8312	DMS7312	DMS6312	DMS5312	DMS4312	DMS3312	 CM7, Direct ,3m
New generation (Masim)	DMS9312A	DMS8312A	DMS7312A	DMS6312A	DMS5312A	DMS4312A	DMS3312A	 CM7, Direct ,3m
<b>Jeteem</b>								
Bd9000	DMS9349	DMS8349	DMS7349	DMS6349	DMS5349	DMS4349	DMS3349	 CP6R, Direct ,3m
<b>Kontron</b>								
7840, 7845, Minimom 7137, Micromon,Supermon, Kolormon, Trakmon	DMS9313 0608010	DMS8313	DMS7313	DMS6313	DMS5313	DMS4313	DMS3313	 CN20, Direct ,3m
7 pins 80 degree	DMS9313A	DMS8313A	DMS7313A	DMS6313A	DMS5313A	DMS4313A	DMS3313A	 CP7T, Direct ,3m
<b>Konica Minolta</b>								
Portable Pulsox-300, Pulsox-300i	DMS9177	DMS8177	DMS7177	DMS6177	DMS5177	DMS4177	DMS3177	 CN100, Direct ,3m
<b>L&amp;T Medical</b>								
Stellar P, Star 50, Galaxy, Planet, Planet 50, Cosmic	DMS9376	DMS8376	DMS7376	DMS6376	DMS5376	DMS4376	DMS3376	 CN79, Direct ,3m
	DMS9376A	DMS8376A	DMS7376A	DMS6376A	DMS5376A	DMS4376A	DMS3376A	
Stellar P, Star 50, Galaxy, Planet, Planet 50, Cosmic	DMS9376B	DMS8376B	DMS7376B	DMS6376B	DMS5376B	DMS4376B	DMS3376B	 CP8R, Direct ,3m

Reusable SpO2 Sensors

Monitor	Adult clip	Adult soft	Pedia. Clip	Pedia. Soft	Neonate wrap	Multi-Y	Adult Ear clip	Connector/Adapter/Len.
<b>Masim</b>								
GE PRO1000,V300, NVIVO:T8l, Schiller,Datascope, ZOLL R Series, MDL Lifepak20 (Masim tech)	DMS9315 1396/1969	DMS8315	DMS7315	DMS6315	DMS5315	DMS4315	DMS3315	 CN22, Direct ,3m
(Masim tech)	DMS9115A	DMS8115A	DMS7115A	DMS6115A	DMS5115A	DMS4115A	DMS3115A	 CN23, See Masim cable, 1m
DB9PIN (Masim Oximax)	DMS9115F 1863/LNCS DCI	DMS8115F	DMS7115F	DMS6115F	DMS5115F	DMS4115F	DMS3115F	
(Masim tech)	DMS9115B LNOP DCI	DMS8115B	DMS7115B LNOP DC-IP	DMS6115B	DMS5115B	DMS4115B 1544	DMS3115B LNOP TC-I	 CN24, See Masim cable, 1m
(Masim Oximax)	DMS9115K	DMS8115K	DMS7115K	DMS6115K	DMS5115K	DMS4115K	DMS3115K	
(Masim Oximax)	DMS9115H	DMS8115H	DMS7115H	DMS6115H	DMS5115H	DMS4115H	DMS3115H	 CN86, See Masim cable, 1m
Radical-7,Rad-87,Rad-57, Pronto, RED-8,ZOLL X Series (Masim Oximax)	DMS9315G	DMS8315G	DMS7315G	DMS6315G	DMS5315G	DMS4315G	DMS3315G	 CN80, Direct ,3m
Corplus , Philips Efficia CM100 , Zoll (Masim Oximax) (N)	DMS9315Q	DMS8315Q	DMS7315Q	DMS6315Q	DMS5315Q	DMS4315Q	DMS3315Q	 CN441, Direct ,3m
RD rainbow,CM 100/120/150 Corplus ,Philips Efficia CM100,Zoll(Masim Oximax) (N)	DMS9315Q-B	DMS8315Q-B	DMS7315Q-B	DMS6315Q-B	DMS5315Q-B	DMS4315Q-B	DMS3315Q-B	 CN441-B, Direct ,3m
RD SET (Masim Oximax) (N)	DMS9115S RD SET DC-I	DMS8115S RD SET DB-I	DMS7115S RD SET DCI-P	DMS6115S	DMS5115S	DMS4115S RD SET Y-I Multisite Sensor	DMS3115S RD SET TC-I	 CN443, Direct ,1m
Rad-97 (N)	DMS9315T	DMS8315T	DMS7315T	DMS6315T	DMS5315T	DMS4315T	DMS3315T	 CN457, Direct ,3m
<b>Medlab</b>								
NANOX2, NANOXeco	DMS9303A	DMS8303A	DMS7303A	DMS6303A	DMS5303A	DMS4303A	DMS3303A	 CN81, Direct ,3m
PEARL10, NANOX10, Nanox10C	DMS9319	DMS8319	DMS7319	DMS6319	DMS5319	DMS4319	DMS3319	 CN30, Direct ,3m
CAPNOX, VITRO, élan, POX 10, P-OX100L	DMS9303	DMS8303	DMS7303	DMS6303	DMS5303	DMS4303	DMS3303	 CP7, Direct ,3m

Reusable SpO2 Sensors

Monitor	Adult clip	Adult soft	Pedia. Clip	Pedia. Soft	Neonate wrap	Multi-Y	Adult Ear clip	Connector/Adapter/Len.
<b>Medsorce</b>								
BPM-200	DMS9319	DMS8319	DMS7319	DMS6319	DMS5319	DMS4319	DMS3319	 CN30, Direct ,3m
<b>MEK</b>								
MP 100,110,400,1000	DMS9316	DMS8316	DMS7316	DMS6316	DMS5316	DMS4316	DMS3316	 CN30B, Direct ,3m
	DMS9116	DMS8116	DMS7116	DMS6116	DMS5116	DMS4116	DMS3116	 CN30B, Direct ,1m
Normal 8 pins	DMS9316A	DMS8316A	DMS7316A	DMS6316A	DMS5316A	DMS4316A	DMS3316A	 CN65, Direct ,3m
Mp1000	DMS9316B	DMS8316B	DMS7316B	DMS6316B	DMS5316B	DMS4316B	DMS3316B	
PAMO II	DMS9316C	DMS8316C	DMS7316C	DMS6316C	DMS5316C	DMS4316C	DMS3316C	 CN66, Direct ,3m
<b>Mennen</b>								
Horizon XL, Mercury 1200,MR 1300, Envoy	DMS9317	DMS8317	DMS7317	DMS6317	DMS5317	DMS4317	DMS3317	 CM10, Direct ,3m
Normal 8 pins	DMS9317A	DMS8317A	DMS7317A	DMS6317A	DMS5317A	DMS4317A	DMS3317A	 CN25, Direct ,3m
(Masim tech)	DMS9317B	DMS8317B	DMS7317B	DMS6317B	DMS5317B	DMS4317B	DMS3317B	 CN26, Direct ,3m
(Oximax tech)	DMS9317D	DMS8317D	DMS7317D	DMS6317D	DMS5317D	DMS4317D	DMS3317D	
<b>Mindray</b>								
12 Male pins	DMS9318	DMS8318	DMS7318	DMS6318	DMS5318	DMS4318	DMS3318	 CN94, Direct ,3m
12 Male pins (Masim tech)	DMS9318A	DMS8318A	DMS7318A	DMS6318A	DMS5318A	DMS4318A	DMS3318A	
503,507,508,509 (Criticare tech)	DMS9318B	DMS8318B	DMS7318B	DMS6318B	DMS5318B	DMS4318B	DMS3318B	 CM5, Direct ,3m
Normal 5 pins	DMS9318C	DMS8318C	DMS7318C	DMS6318C	DMS5318C	DMS4318C	DMS3318C	 CP5, Direct ,3m

Reusable SpO2 Sensors

Monitor	Adult clip	Adult soft	Pedia. Clip	Pedia. Soft	Neonate wrap	Multi-Y	Adult Ear clip	Connector/Adapter/Len.
<b>Mindray</b>								
Normal 6 pins	DMS9318D	DMS8318D	DMS7318D	DMS6318D	DMS5318D	DMS4318D	DMS3318D	 CP6, Direct ,3m
Mindray :PM7000, Pm8000, PM9000,MEC-1000,MEC-2000 Creative: UP-6000, Pc9000 Edan:M9,M99A,M8, M8A,M8B,M3,M3B	DMS9318E	DMS8318E	DMS7318E	DMS6318E	DMS5318E	DMS4318E	DMS3318E	 CP6R, Direct ,3m
6 male pins, 40 degree <b>(Masim tech)</b>	DMS9318F	DMS8318F	DMS7318F	DMS6318F	DMS5318F	DMS4318F	DMS3318F	 CP6S, Direct ,3m
6 male pins, 60 degree <b>(Masim tech)</b>	DMS9318G	DMS8318G	DMS7318G	DMS6318G	DMS5318G	DMS4318G	DMS3318G	 CP6T, Direct ,3m
6 male pins, 80 degree <b>(Nellcor Oximax)</b>	DMS9318L	DMS8318L	DMS7318L	DMS6318L	DMS5318L	DMS4318L	DMS3318L	 CN27, Direct ,3m
PM-6800,BeneView T5 OEM No.0010-30-42737, VS600/VS900,iMEC/iMEC 12, Penlon SP M8,EDAN X12	DMS9318H	DMS8318H	DMS7318H	DMS6318H	DMS5318H	DMS4318H	DMS3318H	 CN28, Direct ,3m
T5/T8 <b>(Masim tech)</b>	DMS9318J	DMS8318J	DMS7318J	DMS6318J	DMS5318J	DMS4318J	DMS3318J	 CN91, Direct ,3m
<b>(Masim Oximax),DPM</b>	DMS9318P	DMS8318P	DMS7318P	DMS6318P	DMS5318P	DMS4318P	DMS3318P	 CN29, Direct ,3m
T5/T8 <b>(Nellcor Oximax)</b>	DMS9318K	DMS8318K	DMS7318K	DMS6318K	DMS5318K	DMS4318K	DMS3318K	 CN07, Direct ,3m
<b>Mindray-Datascope</b>								
Accusat, Accutor 3/4/ SAT/PLUSPassport <b>(Datascope Tech)</b>	DMS9308	DMS8308	DMS7308	DMS6308	DMS5308	DMS4308	DMS3308	 CN08, P0208, 1m
Accusat, Accutor 3/4/ SAT/PLUSPassport <b>(Datascope Tech)</b>	DMS9108A	DMS8108A	DMS7108A	DMS6108A	DMS5108A	DMS4108A	DMS3108A	 CN84, Direct ,3m
Passport II/IILT, XG, Accutor Plus, Expert <b>(Nellcor)</b> <b>(Trio/Spectrum w/opt. Nellcor)</b>	DMS9319B	DMS8319B	DMS7319B	DMS6319B	DMS5319B	DMS4319B	DMS3319B	 CN22, Direct ,3m
Passport II/IILT, XG, Accutor Plus,Trio, Spectrum <b>(Masim)</b>	DMS9315	DMS8315	DMS7315	DMS6315	DMS5315	DMS4315	DMS3315	 CN85, Direct ,3m

Reusable SpO2 Sensors

Monitor	Adult clip	Adult soft	Pedia. Clip	Pedia. Soft	Neonate wrap	Multi-Y	Adult Ear clip	Connector/Adapter/Len.
<b>Mindray-Datascope</b>								
Passport II/IILT, XG, Accutor Plus,Trio, Spectrum <b>(Masim)</b>	DMS9315B	DMS8315B	DMS7315B	DMS6315B	DMS5315B	DMS4315B	DMS3315B	 CN24, Direct ,3m
<b>M&amp;B</b>								
CD531, CD2000,MB526T	DMS9351	DMS8351	DMS7351	DMS6351	DMS5351	DMS4351	DMS3351	 CN70, Direct ,3m
<b>Nellcor</b>								
NPB 40, NPB20, NPB75, NPB 180, NPB185, NPB-190, NPB-195, NPB 3900, NPB 4000	DMS9319	DMS8319	DMS7319	DMS6319	DMS5319	DMS4319	DMS3319	 CN30, Direct ,3m
	DMS9119 512F-30-28263	DMS8119 512E-30-90391	DMS7119 512H-30-79061	DMS6119	DMS5119 518B-30-72107	DMS4119	DMS3119	 CN30, See cable, 1m
N550,N560,N595,N600, N751, N5500 DASH 2500, SOLAR 8000, 8000M,N65 Except CARESCAPE V100	DMS9319A	DMS8319A	DMS7319A	DMS6319A	DMS5319A	DMS4319A	DMS3319A	 CN31, Direct ,3m
	DMS9119A DS-100A	DMS8119A	DMS7119A	DMS6119A	DMS5119A	DMS4119A D-YS	DMS3119A	 CN31, See Nellcor, 1m
NPB3000-3100, Symphony, NPB290,295,395	DMS9119B	DMS8119B	DMS7119B	DMS6119B	DMS5119B	DMS4119B	DMS3119B	 CN22, Direct ,3m
DOC-10 direct sensor Mediana YM-2002 Series	DMS9119C	DMS8119C	DMS7119C	DMS6119C	DMS5119C	DMS4119C	DMS3119C	 CN46, Direct ,3m
<b>Nihon Kohden</b>								
AL-800PA Module; BSM-1101, 2101/2, 3101,7100, 7200; OLV-1100A, Lifescope	DMS9321 TL-101S	DMS8321	DMS7321	DMS6321	DMS5321	DMS4321	DMS3321	 CN32, Direct ,3m
BSM-2301K, BSM-2303K, BSM-2351,BSM-2353, BSM-4103, BSM-4113, VISMO,PVMZ701	DMS9121A TL-201T	DMS8121A	DMS7121A	DMS6121A	DMS5121A TL-220T	DMS4121A	DMS3121A	 CN33, P0221B, 1m
	DMS9321B	DMS8321B	DMS7321B	DMS6321B	DMS5321B	DMS4321B	DMS3321B	 CN34, Direct ,3m
Lifescope OPV 1500 1500k	DMS9321C	DMS8321C	DMS7321C	DMS6321C	DMS5321C	DMS4321C	DMS3321C	 CN84, Direct ,3m
BSM Series <b>(Nellcor Oximax)</b>	DMS9321D	DMS8321D	DMS7321D	DMS6321D	DMS5321D	DMS4321D	DMS3321D	 CN85, Direct ,3m
BSM Series <b>(Masim)</b>	DMS9321E	DMS8321E	DMS7321E	DMS6321E	DMS5321E	DMS4321E	DMS3321E	 CN85, Direct ,3m

Reusable SpO2 Sensors

Monitor	Adult clip	Adult soft	Pedia. Clip	Pedia. Soft	Neonate wrap	Multi-Y	Adult Ear clip	Connector/Adapter/Len.
<b>Nihon Kohden</b>								
SVM-7601	DMS9321F	DMS8321F	DMS7321F	DMS6321F	DMS5321F	DMS4321F	DMS3321F	 CP9R, Direct ,3m
<b>Nonin</b>								
8604	DMS9322	DMS8322	DMS7322	DMS6322	DMS5322	DMS4322	DMS3322	 CM6, Direct ,3m
8500, PamSAT 2500, Avant 9600, Avant 9700, Avant 2120, 9840	DMS9122B 8000AA-1+8500VI	DMS8122B	DMS7122B	DMS6122B	DMS5122B	DMS4122B	DMS3122B	 CN83, Direct ,1m
<b>Novamatrix</b>								
505, 510, 511, 515B, 515C, 520, 7000, 7100, Oxypleth	DMS9323 8776-00	DMS8323	DMS7323	DMS6323	DMS5323	DMS4323	DMS3323	 CN36, Direct ,3m
Dixtal (Brazil)	DMS9323C	DMS8323C	DMS7323C	DMS6323C	DMS5323C	DMS4323C	DMS3323C	 CN36B, Direct ,3m
512/513/710/715,MARS	DMS9123A	DMS8123A	DMS7123A	DMS6123A	DMS5123A	DMS4123A	DMS3123A	 CN37, Direct ,1m
512/513/710/715,MARS	DMS9123B	DMS8123B	DMS7123B	DMS6123B	DMS5123B	DMS4123B	DMS3123B	 CN04, P0223, 1m
<b>Neusoft</b>								
Nova-3D	DMS9343	DMS8343	DMS7343	DMS6343	DMS5343	DMS4343	DMS3343	 CP7, Direct ,3m
General(Digital) Newtech NT3A	DMS9343A	DMS8343A	DMS7343A	DMS6343A	DMS5343A	DMS4343A	DMS3343A	 CN39 Direct ,3m
(Nellcor Oximax)	DMS9343A	DMS8343A	DMS7343A	DMS6343A	DMS5343A	DMS4343A	DMS3343A	 CP8, Direct ,3m
(Digital)	DMS9343C	DMS8343C	DMS7343C	DMS6343C	DMS5343C	DMS4343C	DMS3343C	 CP6R, Direct ,3m
<b>NNC Medical</b>								
(BCI)	DMS9370	DMS8370	DMS7370	DMS6370	DMS5370	DMS4370	DMS3370	 CM6-OB, Direct ,3m
(BCI)	DMS9370A	DMS8370A	DMS7370A	DMS6370A	DMS5370A	DMS4370A	DMS3370A	 CN92, Direct ,3m
Zoncare PM-7000M (Digital)	DMS9370B	DMS8370B	DMS7370B	DMS6370B	DMS5370B	DMS4370B	DMS3370B	

Reusable SpO2 Sensors

Monitor	Adult clip	Adult soft	Pedia. Clip	Pedia. Soft	Neonate wrap	Multi-Y	Adult Ear clip	Connector/Adapter/Len.
<b>Northern Meditech</b>								
Northern Meditech: Gemini, Taurus, Virgo, Aquarius Progetti (Italy): pg S10 aXcent (Germany): GETUS X12	DMS9392	DMS8392	DMS7392	DMS6392	DMS5392	DMS4392	DMS3392	 CN98, Direct ,3m
<b>Nihon Kohden</b>								
NELLCOR OXIMAX	DMS9397	DMS8397	DMS7397	DMS6397	DMS5397	DMS4397	DMS3397	 CN22, Direct ,3m
<b>Osen</b>								
Osen	DMS9359	DMS8359	DMS7359	DMS6359	DMS5359	DMS4359	DMS3359	 CP6R, Direct ,3m
	DMS9359AD	DMS8359AD	DMS7359AD	DMS6359AD	DMS5359AD	DMS4359AD	DMS3359AD	
9000D	DMS9359B	DMS8359B	DMS7359B	DMS6359B	DMS5359B	DMS4359B	DMS3359B	 CP6, Direct ,3m
<b>Pace Tech Generra</b>								
Vitalmax 800 Series	DMS9359C	DMS8359C	DMS7359C	DMS6359C	DMS5359C	DMS4359C	DMS3311C	 CN19, Direct ,3m
<b>Palco</b>								
All	DMS9324	DMS8324	DMS7324	DMS6324	DMS5324	DMS4324	DMS3324	 CN38, Direct ,3m
Palco	DMS9324A	DMS8324A	DMS7324A	DMS6324A	DMS5324A	DMS4324A	DMS3324A	 CN67, Direct ,3m
<b>Philips-HP- Agilent</b>								
78352C, 78354C, 78834C, M1020A, M1025B, M1025A, M1350B, M1722A/B, M1732A/B, M2475B, ViridaA1, A3	DMS9325 M1900B+ M1196T	DMS8325 M1900B+ M1191T	DMS7325	DMS6325 M1900B+ M1192T	DMS5325 M1900B+ M1193T	DMS4325	DMS3325	 CN39, Direct ,3m
MP20, MP30, MP40/50/60/70/90, VM4, VM6, VM8, M3046A	DMS9325A	DMS8325A	DMS7325A	DMS6325A	DMS5325A	DMS4325A	DMS3325A	 CN40, Direct ,3m
	DMS9125A	DMS8125A	DMS7125A	DMS6125A	DMS5125A	DMS4125A	DMS3125A	
Philips-Masim	DMS9325B	DMS8325B	DMS7325B	DMS6325B	DMS5325B	DMS4325B	DMS3325B	 CN40, Direct ,1m
	DMS9325L	DMS8325L	DMS7325L	DMS6325L	DMS5325L	DMS4325L	DMS3325L	
	DMS9325L	DMS8325L	DMS7325L	DMS6325L	DMS5325L	DMS4325L	DMS3325L	 CN41, Direct ,3m

Reusable SpO2 Sensors/

Monitor	Adult clip	Adult soft	Pedia. Clip	Pedia. Soft	Neonate wrap	Multi-Y	Adult Ear clip	Connector/Adapter/Len.
<b>PETAS</b>								
KMA 800	DMS9379B	DMS8379B	DMS7379B	DMS6379B	DMS5379B	DMS4379B	DMS3379B	 CP6R, Direct 3m
<b>Physio-Control</b>								
Lifepak 12	DMS9319	DMS8319	DMS7319	DMS6319	DMS5319	DMS4319	DMS3319	 CN30, Direct 3m
Lifepak 12, 20 (Masim tech)	DMS9315	DMS8315	DMS7315	DMS6315	DMS5315	DMS4315	DMS3315	 CN22, Direct 3m
<b>Perlong</b>								
Perlong PDJ-3000, EDANINS EMF9000B UN-Medical galaxy Adecon DK-8000S SOMO PM-12 Biocare IM12E MILLION ML-1200 KP-900	DMS9364	DMS8364	DMS7364	DMS6364	DMS5364	DMS4364	DMS3364	 CP6, Direct 3m
<b>Primedic</b>								
Primedic	DMS9375	DMS8375	DMS7375	DMS6375	DMS5375	DMS4375	DMS3375	 CP8R, Direct 3m
<b>Schiller</b>								
RGB Omicrom FT, Omicrom PREMIUM	DMS9303	DMS8303	DMS7303	DMS6303	DMS5303	DMS4303	DMS3303	 CP7, Direct 3m
Schiller Truscope (Nellcor oximax)	DMS9326C	DMS8326C	DMS7326C	DMS6326C	DMS5326C	DMS4326C	DMS3326C	
CM8, Schiller, Argus, TM-7 (Nellcor tech)	DMS9326	DMS8326	DMS7326	DMS6326	DMS5326	DMS4326	DMS3326	 CN42, Direct 3m
TM7, Argus LCM (Masim)	DMS9326A	DMS8326A	DMS7326A	DMS6326A	DMS5326A	DMS4326A	DMS3326A	
TM 910	DMS9326B	DMS8326B	DMS7326B	DMS6326B	DMS5326B	DMS4326B	DMS3326B	 CN43, Direct 3m
<b>Spacelabs</b>								
90419, 90501, 90651A-0390351-0/6, 90465, 90466, 90467, 90651A-08, PC Express, 90489, Im77	DMS9327	DMS8327	DMS7327	DMS6327	DMS5327	DMS4327	DMS3327	 CN44, Direct 3m
90496 Ultraview	DMS9327A	DMS8327A	DMS7327A	DMS6327A	DMS5327A	DMS4327A	DMS3327A	 CN45, Direct 3m

Reusable SpO2 Sensors/

Monitor	Adult clip	Adult soft	Pedia. Clip	Pedia. Soft	Neonate wrap	Multi-Y	Adult Ear clip	Connector/Adapter/Len.
<b>Spacelabs</b>								
Mcare300 (Oximax tech)	DMS9327D	DMS8327D	DMS7327D	DMS6327D	DMS5327D	DMS4327D	DMS3327D	 CN72, Direct 3m
<b>Solaris</b>								
(Digital)	DMS9320	DMS8320	DMS7320	DMS6320	DMS5320	DMS4320	DMS3320	 CP5, Direct 3m
Solaris	DMS9320F	DMS8320F	DMS7320F	DMS6320F	DMS5320F	DMS4320F	DMS3320F	
General G2A, G3A, G3C, G3H, G3F (Digital)	DMS9320A	DMS8320A	DMS7320A	DMS6320A	DMS5320A	DMS4320A	DMS3320A	 CP6R, Direct 3m
(Digital)	DMS9120B	DMS8120B	DMS7320B	DMS6120B	DMS5120B	DMS4120B	DMS3120B	 CN30, Direct 1m
Spurt (Digital)	DMS9320C	DMS8320C	DMS7320C	DMS6320C	DMS5320C	DMS4320C	DMS3320C	 CP7, Direct 3m
Creative UP-7000, PC-5000A (Digital)	DMS9320D	DMS8320D	DMS7320D	DMS6320D	DMS5320D	DMS4320D	DMS3320D	 CP5T, Direct 3m
(Digital)	DMS9320E	DMS8320E	DMS7320E	DMS6320E	DMS5320E	DMS4320E	DMS3320E	 CN40, Direct 3m
<b>Senmind</b>								
KP-9000	DMS9348	DMS8348	DMS7348	DMS6348	DMS5348	DMS4348	DMS3348	 CN40, Direct 3m
(Digital)	DMS9348A	DMS8348A	DMS7348A	DMS6348A	DMS5348A	DMS4348A	DMS3348A	
<b>Tianrong Medical</b>								
TR-900E	DMS9357	DMS8357	DMS7357	DMS6357	DMS5357	DMS4357	DMS3357	 CP5, Direct 3m
Tianrong Medical	DMS9357A	DMS8357A	DMS7357A	DMS6357A	DMS5357A	DMS4357A	DMS3357A	 CP6R, Direct 3m
<b>Yonker</b>								
Waltvick Zafiro	DMS93124	DMS83124	DMS73124	DMS63124	DMS53124	DMS43124	DMS33124	 CN438, Direct 3m



## Reusable SpO2 Sensors
















Monitor	Adult clip	Adult soft	Pedia. Clip	Pedia. Soft	Neonate wrap	Multi-Y	Adult Ear clip	Connector/Adapter/Len.
<b>Welch Allyn Protocol</b>								
Propaq 206-EL, 204-EL, 202-EL, 106-EL	DMS9319	DMS8319	DMS7319	DMS6319	DMS5319	DMS4319	DMS3319	 CN30, Direct ,3m
<b>wisonic</b>								
wisonic	DMS91104	DMS81104	DMS71104	DMS61104	DMS51104	DMS41104	DMS31104	 CN92, Direct ,1m
(Digital)	DMS93104	DMS83104	DMS73104	DMS63104	DMS53104	DMS43104	DMS33104	 CN92, Direct ,3m
<b>Zondan</b>								
Zd120	DMS9345	DMS8345	DMS7345	DMS6345	DMS5345	DMS4345	DMS3345	 CP6, Direct ,3m
<b>UTECH</b>								
UTECH	DMS9127	DMS8127	DMS7127	DMS6127	DMS5127	DMS4127	DMS3127	 CP6, Direct ,3m
<b>Konsung</b>								
Aurora	DMS9133	DMS8133	DMS7133	DMS6133	DMS5133	DMS4133	DMS3133	 CN39, Direct ,3m
<b>KERNEL</b>								
KERNEL	DMS9364B	DMS8364B	DMS7364B	DMS6364B	DMS5364B	DMS4364B	DMS3364B	 CP6, Direct ,3m
<b>hlmedical</b>								
TE-4000X	DMS91128	DMS81128	DMS71128	DMS61128	DMS51128	DMS41128	DMS31128	 CM5, Direct ,1m
<b>CO+LINN</b>								
CO+LINN	DMS93131	DMS83131	DMS73131	DMS63131	DMS53131	DMS43131	DMS33131	 CN93, Direct ,3m



Adapter Cables

SpO2 Adapter Cables

















Product types	Photo	P/N Example	P/N Guide
SpO2 Adapter Cables		DMS0221B	P: SpO2 0: SpO2 Adapter Cables 2: 2.2m (1: 1.0m; 3: 3.0m) 21B: NIHON KOHDEN
Warranty: one-year warranty for reusable sensors and cables Package: 1pcs/bag			

Monitor	Medke P/N	Use with sensor	Monitor side conn. P/N	Sensor side conn. P/N
<b>Artema-S&amp;W</b>	DMS0201	Nellcor non-oximax	 CN01	 CN49
	DMS0201A	Ohmeda	 CN02	
<b>Baxter</b>	DMS0202	Nellcor non-oximax/BCI	 CN03	 CN49
<b>Biolight</b>	DMS0205	Nellcor non-oximax	 CP5	 CN49
	DMS0205A	Digital		
	DMS0205B	Nellcor Oximax	 CP7	
	DMS0205D	Digital	 CN61	
	DMS0205E	BLT Db9 digital	 CN62	
	DMS0205F	Nellcor Oximax	 CN61	
	DMS0205J	BLT new oximax digital	 CN82	
	DMS0205H	Digital		
	DMS0205L	BLT new oximax analog	 CN82B	
<b>Bionet</b>	DMS0206	Nellcor non-oximax	 CP7	 CN49

SpO2 Adapter Cables/血氧转接线

Monitor	Medke P/N	Use with sensor	Monitor side conn. P/N	Sensor side conn. P/N
BCI-Smith	DMS0203	BCI	 CP7	 CN50
	DMS0203A 3311	BCI	 CN81	 CN50
Biosys	DMS0204	Nellcor non-oximax/BCI	 CN05	 CN49
CSI	DMS0207 518LD	Nellcor non-oximax/CSI	 CM5	 CN49
	DMS0207A 518DD	Nellcor non-oximax/CSI	 CN06	 CN50
GE-corometrics	DMS0210 4033CAX	Nellcor non-oximax	 CN11	 CN49
Drager	DMS0209 5720084	Nellcor non-oximax	 CM8	 CN49
	DMS0209A 4533840	Nellcor non-oximax	 CN09	
	DMS02098 3375834 3368433	Nellcor non-oximax	 CN10	
	DMS0209DM MS17522	MASIM LNCS		
	DMS0209E	Nellcor oximax	 CN87	
	DMS0209F	Masim new tech		
	DMS0209C MS17041	MASIM LNOP	 CN87	









SpO2 Adapter Cables/

Monitor	Medke P/N	Use with sensor	Monitor side conn. P/N	Sensor side conn. P/N
GE-Trusignal	DMS0210KT TS-G3	GE trusignal	 CN15	 CN57
	DMS0210MT TS-M3	GE trusignal	 CN18	
	DMS021QS TS-H3	GE trusignal	 CN16	 CN56
	DMS0210PS TS-N3	GE trusignal	 CN12	
GE-Datex Ohmeda	DMS0210B OXY-C3	Datex sensor	 CN12	 CN49
	DMS0210BP OXY-SL3	OXY-TIP	 CN12	 CN54
GE- Hellige	DMS0210D 30344358	Nellcor non-oximax	 CN13	 CN49
GE-Marquette	DMS0210E 2006644-001 E9004GE	Nellcor non-oximax	 CN14	 CN49
	DMS0210F 2021406-001	Nellcor oximax		
	DMS0210G Masim tech	Nellcor non-oximax		
	DMS0210GM 2027263-002 (REF 2016)	Masim LNCS	 CN15	
	DMS0210WM	Masim New LNCS		
	DMS0210S PC-12-GE	Masim LNOP	 CN15	

SpO2 Adapter Cables

Monitor	Medke P/N	Use with sensor	Monitor side conn. P/N	Sensor side conn. P/N
GE-Ohmeda	DMS0210H	Ohmeda	 CN16	 CN49
	DMS0210K OXY-ES3	OXY-TIP	 CN15	 CN54
	DMS0210M OXY-MC3	OXY-TIP	 CN18	
	DMS0210Q OXY-OL3	OXY-TIP	 CN16	
Goldway, Pace Tech, Generra	DMS0211	Nellcor non-oximax/BCI	 CM5	 CN49
	DMS0211A	Nellcor non-oximax/BCI	 CP5	
	DMS0211B	Nellcor non-oximax/BCI	 CP7	
	DMS0211C	Nellcor non-oximax/BCI	 CN19	
Invivo	DMS0212	Nellcor non-oximax	 CM7	 CN49
	DMS0212A Masim tech	Nellcor non-oximax		
	DMS0212AM Masim tech	Masim LNCS		
Kontron	DMS0213	Nellcor non-oximax/BCI	 CN20	 CN49
	DMS0213A	Nellcor non-oximax	 CP7T	

SpO2 Adapter Cables

Monitor	Medke P/N	Use with sensor	Monitor side conn. P/N	Sensor side conn. P/N
Mindray	DMS0218 0515-30-11221	Nellcor non-oximax/BCI	 CN39	 CN49
	DMS0218A Masim tech	Nellcor non-oximax/BCI		
	DMS0218AM 0010-30-12452	Masim LNCS		
	DMS0218B 512A-30-06074	Nellcor non-oximax/BCI	 CM5	
	DMS0218C	Nellcor non-oximax/BCI	 CP5	
	DMS0218D	Nellcor non-oximax/BCI	 CP6	
	DMS0218E 0010-20-42594	Nellcor non-oximax/BCI	 CP6R	
	DMS0218F Masim tech	Nellcor non-oximax/BCI		
	DMS0218FM 0010-30-42602	Masim LNCS	 CP6S	
	DMS0218GM 0010-30-42625	Masim LNCS		
	DMS0218H 0010-20-42710	Nellcor non-oximax/BCI		
	DMS0218J Masim tech	Nellcor non-oximax/BCI	 CN28	
DMS0218M 0010-30-42738	Masim LNCS			

SpO2 Adapter Cables

Monitor	Medke P/N	Use with sensor	Monitor side conn. P/N	Sensor side conn. P/N
<b>Mindray</b>	DMS0218K 0010-20-42712 Oximax	Nellcor oximax	 CN29	 CN49
	DMS0218M	Nellcor oximax	 CP6R	
	DMS0218P 115-020768-00	Masim LNCS	 CN91	
	DMS0218L	Nellcor oximax	 CP6T	
	DMS0218N 9200-30-10707	Masim Lnop	 CP6R	
<b>Masim</b>	DMS0215 Masim tech	Nellcor non-oximax/BCI	 CN22	 CN49
	DMS0215M 1814/LNC-10	Masim LNCS	 CN22	 CN49
	DMS0215B Masim 1005/PC08	Masim LNOP	 CN22	 CN55
	DMS0215L	Masim new LNCS	 CN22	 CN49
	DMS0215N 2525 (M-LNC-10)	Masim MLNCS	 CN22	 CN58
	DMS0215G 2059/RED LNC-10	Masim LNCS	 CN80	 CN49
	DMS0215H	Masim MLNCS	 CN80	 CN58
	DMS0215G-S RD rainbow SET Md20	Masim RD SET	 CN80	 CN448

SpO2 Adapter Cables

Monitor	Medke P/N	Use with sensor	Monitor side conn. P/N	Sensor side conn. P/N
<b>Masim</b>	DMS0215Q RED 25 LNC	Masim LNCS	 CN441	 CN49
	DMS0215Q-S	Masim RD SET	 CN441	 CN448
	DMS0215Q-B RED 25 LNC	Masim LNCS	 CN441B	 CN49
	DMS0215Q-BS	Masim RD SET	 CN441B	 CN448
	DMS0215T	Masim RD Set Rad-97	 CN457	 CN448
	DMS0215T-B	Masim LNCS	 CN457	 CN49
	DMS0215T-H	Masim MLNCS	 CN457	 CN58
<b>Mennen</b>	DMS0217	Nellcor non-oximax	 CM10	 CN49
	DMS0217A	Nellcor non-oximax	 CN25	
	DMS0217B	Nellcor non-oximax	 CN26	
	DMS0217BM	Masim LNCS	 CN26	
	DMS0217C	Masim LNOP	 CN26	
<b>Mindray-Datascope</b>	DMS0208 0012-00-0516-02	Datascope	 CN07	 CN49


SpO2 Adapter Cables

Monitor	Medke P/N	Use with sensor	Monitor side conn. P/N	Sensor side conn. P/N
<b>Mindray-Datascope</b>	DMS0215B 0012-00-1099-01	Masim Lnop	 CN22	 CN55
<b>MEK</b>	DMS0216A	DB9 MEK sensor	 CN65	
	DMS0216B	Nellcor non-oximax/BCI	 CN65	 CN49
	DMS0216C	MEK	 CN66	
<b>Nellcor</b>	DMS0219 DEC-8/EC-8/ DEC-4/EC-4	Nellcor oximax	 CN30	 CN49
	DMS0219A	Nellcor non-oximax/BCI	 CN31	 CN51
	DMS0219A Oximax cable	Nellcor non-oximax/BCI	 CN31	 CN51
	DMS0219B SCP-10/MC-10	Nellcor non-oximax	 CN22	 CN49
	DMS0219C DOC-10	Nellcor oximax	 CN46	 CN51
<b>Nihon Kohden</b>	DMS0221	DB9 N.K sensor	 CN32	
	DMS0221B JL-900P	DB9 N.K sensor	 CN34	 CN52
	DMS0221C	DB9 N.K sensor	 CN84	
	DMS0221F	DB9 N.K sensor	 CP9R	



SpO2 Adapter Cables

Monitor	Medke P/N	Use with sensor	Monitor side conn. P/N	Sensor side conn. P/N
<b>Nihon Kohden</b>	DMS0221E	Masim	 CN85	 CN49
	DMS0221D	Nellcor Oximax		
<b>Nonin</b>	DMS0222	Nonin sensor	 CM6	 CN49
	DMS0222A	Nonin sensor	 CN35	
<b>Novamatrix</b>	DMS0223	Novamatrix sensor	 CN36	 CN49
<b>Philips, HP</b>	DMS0225 M1900B	Nellcor non-oximax	 CN39	 CN49
	DMS0225A M1941A	Philips 8 pins sensor	 CN40	 CN53
	DMS0225F 1847/LNOP MP12 Masim	Masim LNOP	 CN41	 CN55
	DMS0225B	Nellcor non-oximax	 CN41	 CN49
	DMS0225BM 2281/LNC MP10 Masim	Masim LNCS sensor		
	DMS0225L-H	Masim MLNCS	 CN41	 CN58
DMS0225C M1940A	Philips 8 pins sensor	 CN39	 CN53	
DMS0225D M1943A/ M1943AL	Nellcor non-oximax	 CN40	 CN49	

SpO2 Adapter Cables

Monitor	Medke P/N	Use with sensor	Monitor side conn. P/N	Sensor side conn. P/N
Philips, HP	DMS0225E M1943NL	Nellcor oximax	 CN40	 CN51
	DMS0225E-S	Masim RD Set	 CN40	 CN448
Solaris-New Tech	DMS0220	DB9 digital sensor	 CP5	 CN49
	DMS0220A	DB9 digital sensor	 CP6R	
	DMS0220D	DB9 digital sensor	 CP5T	
Schiller	DMS0226	BCI sensor	 CN42	 CN49
	DMS0226A	Nellcor Non-Oximax		
	DMS0226AM	Masim LNCS sensor		
	DMS0226B	BCI sensor		
Spacelabs	DMS0227 700-0002-00	Nellcor non-oximax	 CN44	 CN49
	DMS0227A 700-0030-00	Nellcor non-oximax	 CN45	
	DMS0227B	Nellcor non-oximax	 CN69	
	DMS0227BM 700-0906-00	Masim LNCS	 CN69	

SpO2 Adapter Cables

Monitor	Medke P/N	Use with sensor	Monitor side conn. P/N	Sensor side conn. P/N
Spacelabs	DMS0227D	Nellcor oximax	 CN72	 CN49
	DMS0227C	Nellcor oximax	 CN69	 CN55
Contec	DMS02132	DB9 digit sensor	 CN452	 CN49
Comen	DMS0265D	DB9 digit sensor	 CN93	 CN49
	DMS0261M-N	Masim MLNCS	 CN93	 CN58
Co-linn	DMS02131	Nellcor non-oximax	 CN93	 CN49






## ECG Cable & Leadwires

- Fixed One-piece, Trunk cable, Leadwire types available
- Latex-free
- 3 leads, 5 leads, Snap, Pinch available
- AHA and IEC color coding available
- Different types are compatible with most of patient monitors
- Easy to clean and repair




## ECG Cables

Product types	Photo	P/N Example	P/N Guide
Fixed one-piece ECG cable		DMS5101S	G: ECG cables 5: 5 leads (3: 3 leads) 1: AHA (2: IEC) 01: General (10: GE-DATEX; 12: GE-Marquette) S: fixed Snap (P: Fixed pinch)
Trunk cable		DMS52249H	G: ECG cables 5: 5 leads (3: 3 leads) 2: IEC (1: AHA) 24: Philips-HP (10: GE-DATEX) PH: Philips-HP
Leadwire		DMS521PH	G: ECG cables 5: 5 leads (3: 3 leads) 2: IEC (1: AHA) 1: Pinch (2: snap 5: mini pinch) PH: Philips-HP
<p>Warranty: one-year warranty for ECG cables, trunk cables and leadwires Package: 1pcs/bag</p>			



ECG Cables

Monitor	Leads/label	Fixed Pinch	Fixed Snap	Connector
<b>General 6 pins</b>				
Straight connector	3/AHA	DMS3140P	DMS3140S	 CN101A
	3/IEC	DMS3240P	DMS3240S	
	5/AHA	DMS5140P	DMS5140S	
	5/IEC	DMS5240P	DMS5240S	
Angle connector, Lstyle	3/AHA	DMS3101P	DMS3101S	 CN101
	3/IEC	DMS3201P	DMS3201S	
	5/AHA	DMS5101P	DMS5101S	
	5/IEC	DMS5201P	DMS5201S	
<b>General 6 pins</b>				
With (1k resistor) Straight connector	3/AHA	DMS3141P	DMS3141S	 CN101B
	3/IEC	DMS341P	DMS341S	
	5/AHA	DMS514P	DMS514S	
	5/IEC	DMS5241P	DMS5241S	
<b>GE-Critikon</b>				
Series 8700 - Dinamap + Dinamap MPS Select Dinamap MPS Select Portable	3/AHA	DMS3102P	DMS3102S	 CN101A
	3/IEC	DMS3202P	DMS3202S	
	5/AHA	DMS5102P	DMS5102S	
	5/IEC	DMS5202P	DMS5202S	
	3/AHA	DMS31132P	DMS31132S	 CP7R
	3/IEC	DMS32132P	DMS32132S	
	5/AHA	DMS51132P	DMS51132S	
	5/IEC	DMS52132P	DMS52132S	
<b>Artema-S, W</b>				
Series 8000, Athena,Diascope 1-2	3/AHA	DMS3103P	DMS3103S	 CN102
	3/IEC	DMS3203P	DMS3203S	
	5/AHA	DMS5103P	DMS5103S	
	5/IEC	DMS5203P	DMS5203S	
<b>Bionet</b>				
Bm3	3/AHA	DMS3104P	DMS3104S	 CN103
	3/IEC	DMS3204P	DMS3204S	
	5/AHA	DMS5104P	DMS5104S	
	5/IEC	DMS5204P	DMS5204S	
<b>Bionet(146K)</b>				
Bm5	3/AHA	DMS3149P	DMS3149S	 CN132
	3/IEC	DMS3249P	DMS3249S	
	5/AHA	DMS5149P	DMS5149S	
	5/IEC	DMS5249P	DMS5249S	
<b>Biosys</b>				
Petas KMA800	3/AHA	DMS3105P	DMS3105S	 CP6
	3/IEC	DMS3205P	DMS3205S	
	5/AHA	DMS5105P	DMS5105S	
	5/IEC	DMS5205P	DMS5205S	
<b>CO+LINN</b>				
	3/AHA	DMS31131P	DMS31131S	 CN186
	3/IEC	DMS32131P	DMS32131S	
	5/AHA	DMS51131P	DMS51131S	
	5/IEC	DMS52131P	DMS52131S	

ECG Cables

Monitor	Leads/label	Fixed Pinch	Fixed Snap	Connector
<b>Colin</b>				
Mental L connector	3/AHA	DMS3107P	DMS3107S	 CN105
	3/IEC	DMS3207P	DMS3207S	
	5/AHA	DMS5107P	DMS5107S	
	5/IEC	DMS5207P	DMS5207S	
BP 88 - 306 - 508 Plastic connector	3/AHA	DMS3106P	DMS3106S	 CN104
	3/IEC	DMS3206P	DMS3206S	
	5/AHA	DMS5106P	DMS5106S	
	5/IEC	DMS5206P	DMS5206S	
<b>Drager-Siemens</b>				
Sirecust 341-342 R Series 200-400-600-700 -900-1200	3/AHA	DMS3108P	DMS3108S	 CN106
	3/IEC	DMS3208P	DMS3208S	
	5/AHA	DMS5108P	DMS5108S	
	5/IEC	DMS5208P	DMS5208S	
SC5000, SC6000, SC7000, SC8000, Sc9000	3/AHA	DMS3131P	DMS3131S	 CN128
	3/IEC	DMS3231P	DMS3231S	
	5/AHA	DMS5131P	DMS5131S	
	5/IEC	DMS5231P	DMS5231S	
<b>Fukuda Denshi</b>				
DS-7000 DS-7100 DS-7200 DS-7300	3/AHA	DMS3109P	DMS3109S	 CN107
	3/IEC	DMS3209P	DMS3209S	
	5/AHA	DMS5109P	DMS5109S	
	5/IEC	DMS5209P	DMS5209S	
DS-5300W DS-5100E DS-7001	3/AHA	DMS3133P	DMS3133S	 CN123
	3/IEC	DMS3233P	DMS3233S	
	5/AHA	DMS5133P	DMS5133S	
	5/IEC	DMS5233P	DMS5233S	
DYNASCOPE DSL-8001	3/AHA	DMS31144P	DMS31144S	 CN461
	3/IEC	DMS32144P	DMS32144S	
	5/AHA	DMS51144P	DMS51144S	
	5/IEC	DMS52144P	DMS52144S	
<b>GE-Datex Ohmeda</b>				
AS/3 · CS/3 · Cardiocap 1 · Cardiocap 2 · Cardiocap 5 · S/5 · Satellite Plus	3/AHA	DMS3110P	DMS3110S	 CN102
	3/IEC	DMS3210P	DMS3210S	
	5/AHA	DMS5110P	DMS5110S	
	5/IEC	DMS5210P	DMS5210S	
<b>GE-Hellige</b>				
Eagle 1000, Def. SCP 900, Cardioserv, SMC,SMS,SMK,SMU,SMV	3/AHA	DMS3111P	DMS3111S	 CN108
	3/IEC	DMS3211P	DMS3211S	
	5/AHA	DMS5111P	DMS5111S	
	5/IEC	DMS5211P	DMS5211S	
<b>GE-Marquette</b>				
Dash 2000/3000/4000, Dinamap PRO 1000,Eagle Systems, MAC-Lab,Solar 3000/8000/ 9500, Tram 100/200/300/400/ 500/600/ 800/Tram AR,Unity Systems	3/AHA	DMS3112P	DMS3112S	 CN109
	3/IEC	DMS3212P	DMS3212S	
	5/AHA	DMS5112P	DMS5112S	
	5/IEC	DMS5212P	DMS5212S	











ECG Cables

Monitor	Leads/label	Fixed Pinch	Fixed Snap	Connector
<b>Huntleigh Healthcare</b>				
Huntleigh Healthcare	3/AHA	DMS3142P	DMS3142S	 CP8
	3/IEC	DMS3242P	DMS3242S	
	5/AHA	DMS5142P	DMS5142S	
	5/IEC	DMS5242P	DMS5242S	
<b>Kontron</b>				
Micromon, Minimon,Vitalmon, Supermon,Kolormon, Larsen&Toubro Star50	3/AHA	DMS3113P	DMS3113S	 CN110
	3/IEC	DMS3213P	DMS3213S	
	5/AHA	DMS5113P	DMS5113S	
	5/IEC	DMS5213P	DMS5213S	
105, 107, 108, 125, 128, VSM510 Old generation	3/AHA	DMS3114P	DMS3114S	 CN36
	3/IEC	DMS3214P	DMS3214S	
	5/AHA	DMS5114P	DMS5114S	
	5/IEC	DMS5214P	DMS5214S	
<b>Medtronic-Physio Control</b>				
Lifepak 15, 20	3/AHA	DMS3115P	DMS3115S	 CN111C
	3/IEC	DMS3215P	DMS3215S	
	5/AHA	DMS5115P	DMS5115S	
	5/IEC	DMS5215P	DMS5215S	
	10/AHA	DMS3115P	DMS3115S	
<b>Mediana</b>				
Mediana	3/AHA	DMS3129P	DMS3129S	 CP8
	3/IEC	DMS3229P	DMS3229S	
	5/AHA	DMS5129P	DMS5129S	
	5/IEC	DMS5229P	DMS5229S	
<b>MEK</b>				
MP500,MP600,MP1000	3/AHA	DMS3119P	DMS3119S	 CN30
	3/IEC	DMS3219P	DMS3219S	
	5/AHA	DMS5119P	DMS5119S	
	5/IEC	DMS5219P	DMS5219S	
6 pins	3/AHA	DMS3120P	DMS3120S	 CN101B
	3/IEC	DMS3220P	DMS3220S	
	5/AHA	DMS5120P	DMS5120S	
	5/IEC	DMS5220P	DMS5220S	
MEK	3/AHA	DMS3152P	DMS3152S	 CN133
	3/IEC	DMS3252P	DMS3252S	
	5/AHA	DMS5152P	DMS5152S	
	5/IEC	DMS5252P	DMS5252S	
<b>Mennen</b>				
Series 9000-1100-2000	3/AHA	DMS3116P	DMS3116S	 CN112
	3/IEC	DMS3216P	DMS3216S	
	5/AHA	DMS5116P	DMS5116S	
	5/IEC	DMS5216P	DMS5216S	
Envoy, Envoy 1000, Vitalogik 6000	3/AHA	DMS3117P	DMS3117S	 CN113
	3/IEC	DMS3217P	DMS3217S	
	5/AHA	DMS5117P	DMS5117S	
	5/IEC	DMS5217P	DMS5217S	

ECG Cables

Monitor	Leads/label	Fixed Pinch	Fixed Snap	Connector
<b>Mindray</b>				
6 pins (with resistor) PM-7000, PM-8000series, PM-9000series, new PM-6000	3/AHA	DMS3141P	DMS3141S	 CN101A
	3/IEC	DMS3241P	DMS3241S	
	5/AHA	DMS5141P	DMS5141S	
	5/IEC	DMS5241P	DMS5241S	
6 pins (without resistor) MEC1000, MEC1200 Series 9000-1100-2000	3/AHA	DMS3140P	DMS3140S	 CN101A
	3/IEC	DMS3240P	DMS3240S	
	5/AHA	DMS5140P	DMS5140S	
	5/IEC	DMS5240P	DMS5240S	
12 pins (without resistor) PM-5000	3/AHA	DMS3118P	DMS3118S	 CN127
	3/IEC	DMS3218P	DMS3218S	
	5/AHA	DMS5118P	DMS5118S	
	5/IEC	DMS5218P	DMS5218S	
12 pins (with resistor) BENEVIEW T5/T8	3/AHA	DMS3143P	DMS3143S	 CN127
	3/IEC	DMS3243P	DMS3243S	
	5/AHA	DMS5143P	DMS5143S	
	5/IEC	DMS5243P	DMS5243S	
<b>Mindray-Datascope</b>				
V12/21	3/AHA	DMS3145P	DMS3145S	 CN130
	3/IEC	DMS3245P	DMS3245S	
	5/AHA	DMS5145P	DMS5145S	
	5/IEC	DMS5245P	DMS5245S	
<b>Nihon Kohden</b>				
Def. Cardiolife Lifescope series 6-8-9-11 6102-6205-6120	3/AHA	DMS3121P	DMS3121S	 CN114
	3/IEC	DMS3221P	DMS3221S	
	5/AHA	DMS5121P	DMS5121S	
	5/IEC	DMS5221P	DMS5221S	
BSM-2300 lifescope I,BSM-4100 lifescope P,BSM-5100 lifescope A,BSM-9510lifescope M,BSM-9800 lifescope S,VPMZ701	3/AHA	DMS3122P	DMS3122S	 CN115
	3/IEC	DMS3222P	DMS3222S	
	5/AHA	DMS5122P	DMS5122S	
	5/IEC	DMS5222P	DMS5222S	
TEC-5521	3/AHA	DMS3130P	DMS3130S	 CN119
	3/IEC	DMS3218P	DMS3218S	
	5/AHA	DMS5130P	DMS5130S	
	5/IEC	DMS5230P	DMS5230S	
OPV-1500	3/AHA	DMS3134P	DMS3134S	 CN126
	3/IEC	DMS3234P	DMS3234S	
	5/AHA	DMS5134P	DMS5134S	
	5/IEC	DMS5234P	DMS5234S	
<b>Spacelabs</b>				
Ultraview New generation	3/AHA	DMS3126P	DMS3126S	 CN118
	3/IEC	DMS3226P	DMS3226S	
	5/AHA	DMS5126P	DMS5126S	
	5/IEC	DMS5226P	DMS5226S	
<b>Schiller</b>				
Schiller,Bruker,ODM	3/AHA	DMS3125P	DMS3125S	 CN110
	3/IEC	DMS3225P	DMS3225S	
	5/AHA	DMS5125P	DMS5125S	
	5/IEC	DMS5225P	DMS5225S	










ECG Cables

Monitor	Leads/label	Fixed Pinch	Fixed Snap	Connector
<b>Philips-HP</b>				
43100A,43110A,43110M/MC,43120A,43200A,77000AC,77020A/AC/AR,77025A,77030A,77035A,77040A,78208A,78210A,78211A,78330A,78331A,78332A,78333A,78341A,78342A,78345A,78346A,78351T,78352A,78353A,78353B,78354A,78532A,78532B,78534A/B/C,78619A,78620A,78660A/B,78661A/B,78670A,78671A,M1722A/B,M1723A/B,M2400A,M2406A,M2745B,M3500B,M5500B	3/AHA	DMS3123P	DMS3123S	 CN116
	3/IEC	DMS3223P	DMS3223S	
	5/AHA	DMS5123P	DMS5123S	
	5/IEC	DMS5223P	DMS5223S	
	3/AHA	DMS3124P	DMS3124S	 CN117
	3/IEC	DMS3224P	DMS3224S	
	5/AHA	DMS5124P	DMS5124S	
	5/IEC	DMS5224P	DMS5224S	
<b>Biolight</b>				
M7000,M9000,M9500	3/AHA	DMS3135P	DMS3135S	 CN110
	3/IEC	DMS3235P	DMS3235S	
	5/AHA	DMS5135P	DMS5135S	
	5/IEC	DMS5235P	DMS5235S	
A8	3/AHA	DMS3147P	DMS3147S	 CN131
	3/IEC	DMS3247P	DMS3247S	
	5/AHA	DMS5147P	DMS5147S	
	5/IEC	DMS5247P	DMS5247S	
M800	3/AHA	DMS3117P	DMS3117S	 CN314B
	3/IEC	G32117P	DMS32117S	
Q5	3/AHA	DMS31133P	DMS31133S	 CN131
	3/IEC	DMS32133P	DMS32133S	
	5/AHA	DMS51133P	DMS51133S	
	5/IEC	DMS52133P	DMS52133S	
<b>Drager</b>				
10Pins	3/AHA	DMS3116P	DMS3116S	 CN437
	3/IEC	DMS3216P	DMS3216S	
	5/AHA	DMS5116P	DMS5116S	
	5/IEC	DMS5216P	DMS5216S	
<b>Comen</b>				
6Pins	3/AHA	DMS3132P	DMS3132S	 CN101A
	3/IEC	DMS3232P	DMS3232S	
	5/AHA	DMS5132P	DMS5132S	
	5/IEC	DMS5232P	DMS5232S	
12Pins	3/AHA	DMS3155P	DMS3155S	 CN186
	3/IEC	DMS3255P	DMS3255S	
	5/AHA	DMS5155P	DMS5155S	
	5/IEC	DMS5255P	DMS5255S	
<b>SAADAT</b>				
NEPTUNE, HEYER	3/AHA	DMS31121P	DMS31121S	 CN360
	3/IEC	DMS32121P	DMS32121S	
	5/AHA	DMS51121P	DMS51121S	
	5/IEC	DMS52121P	DMS52121S	

ECG Cables

Monitor	Leads/label	Fixed Pinch	Fixed Snap	Connector
<b>Zoll</b>				
10Pins	3/AHA	DMS31119S	DMS31119S	 CN431
	3/IEC	DMS32119S	DMS32119S	
	5/AHA	DMS51119P	DMS51119S	
	5/IEC	DMS52119S	DMS52119S	
<b>Yonker</b>				
9Pins	3/AHA	DMS31118P	DMS31118S	 CN439
	3/IEC	DMS32118S	DMS32118S	
	5/AHA	DMS51118P	DMS51118S	
	5/IEC	DMS52118S	DMS52118S	
<b>PETAS</b>				
KMA 160R,165R,260R,K-250,300	4/AHA	DMS4148P	DMS4148S	 CP6-5F
	4/IEC	DMS4248P	DMS4248S	
	5/AHA	DMS5148P	DMS5148S	
	5/IEC	DMS5248P	DMS5248S	
<b>Creative</b>				
PC9000D	3/AHA	DMS3127P	DMS3127S	 CP6
	3/IEC	DMS3227P	DMS3227S	
	5/AHA	DMS5127P	DMS5127S	
	5/IEC	DMS5227P	DMS5227S	
<b>M&amp;B</b>				
Cd2000	3/AHA	DMS3128P	DMS3128S	 CN121
	3/IEC	DMS3228P	DMS3228S	
	5/AHA	DMS5128P	DMS5128S	
	5/IEC	DMS5228P	DMS5228S	
<b>Infinium</b>				
Infinium	3/AHA	DMS3153P	DMS3153S	 CN185
	3/IEC	DMS3253P	DMS3253S	
	5/AHA	DMS5153P	DMS5153S	
	5/IEC	DMS5253P	DMS5253S	
<b>JincoMed</b>				
11Pins	3/AHA	DMS3199P	DMS3199S	
	3/IEC	DMS3299P	DMS3299S	
	10/AHA	DMS1199P	DMS1199S	
	10/IEC	DMS1299P	DMS1299S	
<b>ARROW</b>				
AUTOCAT 2	3/AHA	DMS31115P	DMS31115S	 CN110
	3/IEC	DMS32115P	DMS32115S	
	5/AHA	DMS51115P	DMS51115S	
	5/IEC	DMS52115P	DMS52115S	
<b>Tian Rong</b>				
6Pins	3/AHA	DMS3154P	DMS3154S	 CP6-5F
	3/IEC	DMS3254P	DMS3254S	
	5/AHA	DMS5154P	DMS5154S	
	5/IEC	DMS5254P	DMS5254S	

## ECG Cables

Monitor	Leads/label	Fixed Pinch	Fixed Snap	Connector
<b>Primedic</b>				
DM10,DM30	4/AHA	DMS4158P	DMS4158S	 CN425
	4/IEC	DMS4258P	DMS4258S	
DefiMonitor EVO IV Mediana D500	3/AHA	DMS3159P	DMS3159S	 CN141
	3/IEC	DMS3259P	DMS3259S	
	5/AHA	DMS5159P	DMS5159S	
	5/IEC	DMS5259P	DMS5259S	
	10/AHA	DMS1159P	DMS1159S	
	10/IEC	DMS1259P	DMS1259S	
Primedic Xd30	3/AHA	DMS31104P	DMS31104S	
	3/IEC	DMS32104P	DMS32104S	
<b>Contec</b>				
CSM9000®	3/AHA	DMS31135P	DMS31135S	 CN453
	3/IEC	DMS32135P	DMS32135S	
	5/AHA	DMS51135P	DMS51135S	
	5/IEC	DMS52135P	DMS52135S	
<b>Bistos</b>				
BT-780®	3/AHA	DMS31124P	DMS31124S	 CN264
	3/IEC	DMS32124P	DMS32124S	
	5/AHA	DMS51124P	DMS51124S	
	5/IEC	DMS52124P	DMS52124S	
<b>NCC /Wisonic</b>				
6Pins	3/AHA	DMS31136P	DMS31136S	 CN138
	3/IEC	DMS32136P	DMS32136S	
	5/AHA	DMS51136P	DMS51136S	
	5/IEC	DMS52136P	DMS52136S	
<b>General</b>				
See Philips	3/AHA	DMS3124P	DMS3124S	 CN117
	3/IEC	DMS3224P	DMS3224S	
	5/AHA	DMS5124P	DMS5124S	
	5/IEC	DMS5224P	DMS5224S	
<b>Narcotrend</b>				
9Pins	3/AHA	DMS31126P	DMS31126S	 CP9R
	3/IEC	DMS32126P	DMS32126S	
	5/AHA	DMS51126P	DMS51126S	
	5/IEC	DMS52126P	DMS52126S	
<b>hlmedical</b>				
TE-4000X	3/AHA	DMS31143P	DMS31143S	
	3/IEC	DMS32143P	DMS32143S	
	5/AHA	DMS51143P	DMS51143S	
	5/IEC	DMS52143P	DMS52143S	


















## ECG Trunk Cable & Leadwires

- Fixed One-piece, Trunk cable, Leadwire types available
- Latex-free
- 3 leads, 5 leads, Snap, Pinch available
- AHA and IEC color coding available
- Different types are Compatible with most of patient monitors
- Easy to clean and repair



Product types	Photo	P/N Example	P/N Guide
Trunk cable		DMS5254PH	G: ECG cables 5: 5 leads (3: 3 leads) 2: IEC (1: AHA) 24: Philips-HP (10: GE-DATEX) PH: Philips-HP
Leadwire		DMS521PH	G: ECG cables 5: 5 leads (3: 3 leads) 2: IEC (1: AHA) 1: Pinch (2: snap 5: mini pinch) PH: Philips-HP
Warranty: one-year warranty for ECG cables, trunk cables and leadwires Package: 1pcs/bag			

ECG Trunk Cable
















Monitor	Leads/label	Cable P/N	Connector	Pinch	Snap	Trunk cable Conn.
<b>General 6 pins</b>						
Straight connector	3/AHA	DMS3140DN		DMS311DN	DMS312DN	
	3/IEC	DMS3240DN		DMS321DN	DMS322DN	
	5/AHA	DMS5140DN		DMS511DN	DMS512DN	
	5/IEC	DMS5240DN		DMS521DN	DMS522DN	
	3/AHA	DMS3140LL		DMS311LL	DMS312LL	
	3/IEC	DMS3240LL		DMS321LL	DMS322LL	
	5/AHA	DMS5140LL		DMS511LL	DMS512LL	
	5/IEC	DMS5240LL		DMS521LL	DMS522LL	
Angle connector	3/AHA	DMS3101DN		DMS311DN	DMS312DN	
	3/IEC	DMS3201DN		DMS321DN	DMS322DN	
	5/AHA	DMS5101DN		DMS511DN	DMS512DN	
	5/IEC	DMS5201DN		DMS521DN	DMS522DN	
LL style	3/AHA	G3101LL		DMS311LL	DMS312LL	
	3/IEC	G3201LL		DMS321LL	DMS322LL	
	5/AHA	G5101LL		DMS511LL	DMS512LL	
	5/IEC	G5201LL		DMS521LL	DMS522LL	
<b>General 6 pins</b>						
Normal 6 pins conn. With (1k resistor) Straight connector	3/AHA	DMS3141DN		DMS311DN	DMS312DN	
	3/IEC	DMS3241DN		DMS321DN	DMS322DN	
	5/AHA	DMS5141DN		DMS511DN	DMS512DN	
	5/IEC	DMS5241DN		DMS521DN	DMS522DN	
<b>GE-Critikon</b>						
Series 8700 - Dinamap +Dinamap MPS Select Dinamap MPS Select Portable	3/AHA	DMS3102DX		DMS311DX	DMS312DX	
	3/IEC	DMS3202DX		DMS321DX	DMS322DX	
	5/AHA	DMS5102DX		DMS511DX	DMS512DX	
	5/IEC	DMS5202DX		DMS521DX	DMS522DX	

ECG Trunk Cable

















Monitor	Leads/label	Cable P/N	Connector	Pinch	Snap	Trunk cable Conn.
<b>Colin</b>						
BP 88 - 306 - 508 Plastic connector	3/AHA	DMS3106CL		DMS311CL	DMS312CL	
	3/IEC	DMS3206CL		DMS321CL	DMS322CL	
	5/AHA	DMS5106CL		DMS511CL	DMS512CL	
	5/IEC	DMS5206CL		DMS521CL	DMS522CL	
<b>Drager-Siemens</b>						
Sirecust 341-342 R Series 200-400-600-700-900-1200	3/AHA	DMS3108DR		DMS311DR	DMS312DR	
	3/IEC	DMS3208DR		DMS321DR	DMS322DR	
	5/AHA	DMS5108DR		DMS511DR	DMS512DR	
	5/IEC	DMS5208DR		DMS521DR	DMS522DR	
Drager-Siemens	3/AHA	DMS3116DR		DMS311DR	DMS312DR	
	3/IEC	DMS3216DR		DMS321DR	DMS322DR	
	5/AHA	DMS5116DR		DMS511DR	DMS512DR	
	5/IEC	DMS5216DR		DMS521DR	DMS522DR	
<b>Fukuda Denshi</b>						
DS-7000 DS-7100 DS-7200 DS-7300	3/AHA	DMS3109DR		DMS311DR	DMS312DR	
	3/IEC	DMS3209DR		DMS321DR	DMS322DR	
	5/AHA	DMS5109DR		DMS511DR	DMS512DR	
	5/IEC	DMS5209DR		DMS521DR	DMS522DR	
<b>GE-Datex Ohmeda</b>						
AS/3, CS/3, Cardiocap 1, Cardiocap 2, Cardiocap 5, S/5, Satellite Plus	3/AHA	DMS3110DX 545302 545307		DMS311DX	DMS312DX	
	3/IEC	DMS3210DX 545305 545300 545304		DMS321DX	DMS322DX	
	5/AHA	DMS5110DX 545303 545308		DMS511DX	DMS512DX	
	5/IEC	DMS5210DX 545301 545306		DMS521DX	DMS522DX	
<b>GE-Hellige</b>						
Eagle 1000, Def SCP 900, Cardioserv, SMC,SMS,SMK, SMU,SMV	3/AHA	DMS3111DX		DMS311DX	DMS312DX	
	3/IEC	DMS3211DX		DMS321DX	DMS322DX	
	5/AHA	DMS5111DX		DMS511DX	DMS512DX	
	5/IEC	DMS5211DX		DMS521DX	DMS522DX	

ECG









ECG Trunk Cable

Monitor	Leads/label	Cable P/N	Connector	Pinch	Snap	Trunk cable Conn.
<b>Mediana</b>						
Mediana	3/AHA	DMS3129DN	 CP8	DMS3129P	DMS3129S	 CN151
	3/IEC	DMS3229DN		DMS3229P	DMS3229S	
<b>GE-Marquette</b>						
Dash 2000/3000/4000, Dinamap PRO 1000, Eagle Systems, MAC-Lab, Solar 3000/8000/9500, Tram 100/200/300/400/ 500/600/800/Tram AR, Unity Systems	3/AHA	DMS5112MQ 412931-001 412931-021	 CN109	DMS311MQ 412682-001 412682-002	DMS312MQ 411203-001 411203-002	 CN167
	3/IEC	DMS5212MQ 2017003-003 412931-002		DMS321MQ 412682-003 412682-004	DMS322MQ 411203-003 411203-004	
	5/AHA	DMS5112MQ		DMS511MQ 412681-011 412681-002	DMS512MQ 412681-003 412681-001	 CN167
	5/IEC	DMS5212MQ		DMS521MQ 412681-011 412681-002	DMS522MQ 412681-003 412681-001	
	6/AHA	DMS6112MQ		 CN136	DMS611MQ	DMS612MQ
	6/IEC	DMS6212MQ			DMS621MQ	DMS622MQ
	10/AHA	DMS1112MQ		 CN251	DMS111MQ	DMS112MQ
	10/IEC	DMS1212MQ			DMS121MQ	DMS122MQ
<b>Mennen</b>						
Series 9000-1100-2000	3/AHA	DMS3116MN	 CN112	DMS311AA	DMS312AA	 CN173
	3/IEC	DMS3216MN		DMS321AA	DMS322AA	
	5/AHA	DMS5116MN		DMS511AA	DMS512AA	 CN174
	5/IEC	DMS5216MN		DMS521AA	DMS522AA	
Envoy, Envoy 1000, Vitalogik 6000	3/AHA	DMS6117MN	 CN113	DMS311AA	DMS312AA	 CN173
	3/IEC	DMS3217MN		DMS321AA	DMS322AA	
	5/AHA	DMS5117MN		DMS511AA	DMS512AA	 CN174
	5/IEC	DMS5217MN		DMS521AA	DMS522AA	
<b>Mindray-Datascope</b>						
V12/21	3/AHA	DMS5145DT	 CN130	DMS311DT	DMS312DT	 CN184
	3/IEC	DMS5245DT		DMS321DT	DMS322DT	
	5/AHA	DMS5145DT		DMS511DT	DMS512DT	
	5/IEC	DMS5245DT		DMS521DT	DMS522DT	









ECG Trunk Cable

Monitor	Leads/label	Cable P/N	Connector	Pinch	Snap	Trunk cable Conn.
<b>Mindray</b>						
6 pins (with resistor) PM-8000series, PM-9000series, new PM-6000	3/AHA	DMS3141MD 0010-30-12257	 CN101A	DMS311MD 0010-30-12263	DMS312MD 0010-30-12267	 CN159
	3/IEC	DMS3241MD 0010-30-12247		DMS321MD 0010-30-12265	DMS322MD 0010-30-12269	
	5/AHA	DMS5141MD 0010-30-12244		DMS511MD 0010-30-12262	DMS512MD 0010-30-12266	 CN160
	5/IEC	DMS5241MD 0010-30-12245		DMS521MD 0010-30-12264	DMS522MD 0010-30-12268	
12 pins (without resistor) PM-5000	3/AHA	DMS3118MD 0010-30-12380	 CN127	DMS311MD 0010-30-12263	DMS312MD 0010-30-12267	 CN159
	3/IEC	DMS3218MD 0010-30-12251		DMS321MD 0010-30-12265	DMS322MD 0010-30-12269	
	5/AHA	DMS5118MD		DMS511MD 0010-30-12262	DMS512MD 0010-30-12266	 CN160
	5/IEC	DMS5218MD		DMS521MD 0010-30-12264	DMS522MD 0010-30-12268	
12 pins with (1K resistor) PM-6000	3/AHA	DMS3124MD 0010-30-12379	 CN117	DMS311MD 0010-30-12263	DMS312MD 0010-30-12267	 CN159
	3/IEC	DMS3224MD 0010-30-12379		DMS321MD 0010-30-12265	DMS322MD 0010-30-12269	
	5/AHA	DMS5124MD 0010-30-12260		DMS511MD 0010-30-12262	DMS512MD 0010-30-12266	 CN160
	5/IEC	DMS5224MD 0010-30-12260		DMS521MD 0010-30-12264	DMS522MD 0010-30-12268	
12 pins with (1K resistor) BENEVIEW T5/T8	3/AHA	DMS3143MD 0010-30-42719	 CN127	DMS3143MD 0010-30-42726	DMS312MD 0010-30-42734	 CN159
	3/IEC	DMS3243MD 0010-30-42719		DMS321MD 0010-30-42725	DMS322MD 0010-30-42733	
	5/AHA	DMS5143MD 0010-30-42719		DMS511MD 0010-30-42727	DMS512MD 0010-30-42735	 CN160
	5/IEC	DMS5243MD 0010-30-42719		DMS521MD 0010-30-42728	DMS522MD 0010-30-42736	
6 pins without MEC1200Series 9000-1100-2000	3/AHA	DMS3140MD 0010-30-12378 0010-30-12242	 CN101A	DMS3110MD	DMS312MD	 CN159
	3/IEC	DMS3240MD 0010-30-12243		DMS321MD	DMS322MD	
<b>Philips-HP</b>						
78352C, 78354C, 78834C / 862472 to 862475 C3 Monitor / 862478, 862479 / 863064, 863065, 863066, 863068 SureSigns VM6, VM8 / M1001A/B, M1002A/B /M3000A, M3001A, M3002A MMS X2 M8102A MP2 Monitor, M8105A MP5 Monitor / M4735A, M3535A/M3536A	3/AHA	DMS3124PH	 CN117	DMS311PH M1671A	DMS312PH M1673A	 CN179
	3/IEC	DMS3224PH M1669A		DMS321PH M1672A	DMS322PH M1674A	
	5/AHA	DMS5124PH		DMS511PH M1968A	DMS512PH M1644A	 CN180
	5/IEC	DMS5224PH M1668A		DMS521PH M1971A	DMS522PH M1645A	
Philips-HP	10/AHA	DMS1124PH M1663A	 CN139	DMS511PH M1976A	DMS512PH M1602A	 CN600
	10/IEC	DMS1224PH M1663A		DMS521PH M1978A	DMS522PH M1644A	











ECG Trunk Cable

Monitor	Leads/label	cable P/N	Connector	Pinch	Snap	Trunk cable Conn.
<b>Philips-HP</b>						
43100A,43110A,43110M/MC, 43120A,43200A,77000AC,77020A/AC/AR,77025A,77030A,77035A,77040A,78208A,78210A,78211A,78330A,78331A,78332A,78333A,78341A,78342A,78345A,78346A,78351T,78352A,78353A,78353B,78354A,78532A,78532B,78534A/B/C,78619A,78620A,78660A/B,78661A/B,78670A,78671A,M1722A/B,M1723A/B,M2400A,M2406A,M2745B,M3500B,M5500B	3/AHA	DMS3123AA M1580A		DMS311AA M1603A	DMS312AA M1605A	
	3/IEC	DMS3223AA M1590A		DMS321AA M1613A	DMS322AA M1615A	
	5/AHA	DMS5123AA M1600A		DMS511AA M1623A	DMS512AA M1625A	
	5/IEC	DMS5223AA M1610A		DMS521AA M1633A	DMS522AA M1635A	
78352C,78354C,78834C,862474 (C3),862475 (C3),862478 (C3),862479 (C3),MP20,MP30,MP40,MP50,MP60,MP70,MP90,M1001A,M1001B,M1002A,M1002B,M1165A,M1166A,M1167A,M1175A,M1177A,M1205A,M1722A,M1722B,M2601A,M2745A,M3000A,M3001A,M3046A,M3500B,M3535A,M3536A,M4375A,M5500B	3/AHA	DMS3124AA M1500A		DMS311AA M1603A	DMS312AA M1605A	
	3/IEC	DMS3224AA M1510A		DMS321AA M1613A	DMS322AA M1615A	
	5/AHA	DMS5124AA M1520A		DMS511AA M1623A	DMS512AA M1625A	
	5/IEC	DMS5224AA M1530A		DMS521AA M1633A	DMS522AA M1635A	
<b>Biolight</b>						
M7000,M9000,M9500	3/AHA	DMS3135DN		DMS311DN	DMS312DN	
	3/IEC	DMS3235DN		DMS321DN	DMS322DN	
	5/AHA	DMS5135DN		DMS511DN	DMS512DN	
	5/IEC	DMS5235DN		DMS521DN	DMS522DN	
<b>Nihon Kohden</b>						
Def. Cardioline Lifescope series 6-8-9-11 6102-6205-6120	3/AHA	DMS3121NH		DMS311NH	DMS312NH	
	3/IEC	DMS3221NH		DMS321NH	DMS322NH	
	5/AHA	DMS5121NH		DMS511NH	DMS512NH	
	5/IEC	DMS5221NH		DMS521NH	DMS522NH	
BSM-2300 lifescope I, BSM-4100 lifescope P, BSM-5100 lifescope A, BSM-9510lifescope M, BSM-9800 lifescope S, VPMZ701	3/AHA	DMS6222NK		DMS311NK	DMS312NK	
	3/IEC	DMS6222NK		DMS321NK	DMS322NK	
	6/AHA	DMS6222NK		DMS611NK	DMS612NK	
	6/IEC	DMS6222NK		DMS621NK	DMS622NK	
TEC-5521	3/AHA	DMS3130NH		DMS311NH	DMS312NH	
	3/IEC	DMS3230NH		DMS321NH	DMS322NH	
OPV-1500	3/AHA	DMS3134NK		DMS311NK	DMS312NK	
	3/IEC	DMS3234NK		DMS321NK	DMS322NK	













ECG Trunk Cable

Monitor	Leads/label	cable P/N	Connector	Pinch	Snap	Trunk cable Conn.	
<b>Creative</b>							
PC9000D	5/AHA	DMS5127DN		DMS511DN	G512DN		
	5/IEC	DMS227DN		DMS21DN	G522DN		
<b>Spacelabs</b>							
Ultraview New generation	3/AHA	DMS5126SL		DMS511SL	DMS512SL		
	3/IEC	DMS226SL		DMS21SL	DMS22SL		
	5/AHA	DMS5126SL		DMS511SL	DMS512SL		
	5/IEC	DMS226SL		DMS21SL	DMS22SL		
<b>SA IRAN (OEM)</b>							
ARAD P10	5/AHA	DMS5141MD-FM		DMS511MD-FM	DMS512MD-FM		
	5/IEC	DMS5241MD-FM		DMS521MD-FM	DMS522MD-FM		
	10/AHA	DMS11109MQ			DMS111MQ	DMS121MQ	
	10/IEC	DMS12109MQ			DMS112MQ	DMS122MQ	

ECG Leadwires

Type Designation	Leads/label	Trunk cable Conn.	Pinch	Snap	Leadwires Conn.
<b>LI type</b>	3/AHA	 CN154	DMS311LL	DMS312LL	 CN156
	3/IEC		DMS321LL	DMS322LL	
	5/AHA	 CN155	DMS511LL	DMS512LL	 CN156
	5/IEC		DMS521LL	DMS522LL	
<b>Din type</b>	3/AHA	 CN151	DMS311DN	DMS312DN	 CN153
	3/IEC		DMS321DN	DMS322DN	
	5/AHA	 CN152	DMS511DN	DMS512DN	 CN153
	5/IEC		DMS521DN	DMS522DN	
<b>Colin</b>	3/AHA	 CN451	DMS311CL	DMS312CL	 CN158B
	3/IEC		DMS321CL	DMS322CL	
	5/AHA		DMS511CL	DMS512CL	 CN158A
	5/IEC		DMS521CL	DMS522CL	

ECG Leadwires

Type Designation	Leads/label	Trunk cable Conn.	Pinch	Snap	Leadwires Conn.
<b>GE-Datex Ohmeda</b>	3/AHA	 CN163	DMS311DX	DMS312DX	 CN165
	3/IEC		DMS321DX	DMS322DX	
	5/AHA	 CN164	DMS511DX	DMS512DX	 CN166
	5/IEC		DMS521DX	DMS522DX	
<b>Philips-HP</b>	3/AHA	 CN173	DMS311DX	DMS312DX	 CN175
	3/IEC		DMS321AA	DMS322AA	
	5/AHA	 CN174	DMS511AA	DMS512AA	 CN176
	5/IEC		DMS521AA	DMS522AA	
	3/AHA	 CN179	DMS311PH	DMS312PH	 CN181
	3/IEC		DMS321PH	DMS322PH	
	5/AHA	 CN180	DMS511PH	DMS512PH	 CN182
	5/IEC		DMS521PH	DMS522PH	

ECG










ECG Leadwires/心电分线

Type Designation 型号类型	Leads/label 导联数	Trunk cable Conn. 分线盒接头	Pinch 夹式	Snap 扣式	Leadwires Conn. 导联线接头
LI type	3/AHA	CN154	G311LL	G312LL	CN156
	3/IEC		G321LL	G322LL	
	5/AHA	CN155	G511LL	G512LL	CN156
	5/IEC		G521LL	G522LL	
Din type	3/AHA	CN151	G311DN	G312DN	CN153
	3/IEC		G321DN	G322DN	
	5/AHA	CN152	G511DN	G512DN	CN153
	5/IEC		G521DN	G522DN	
Colin 科林	3/AHA	CN451	G311CL	G312CL	CN158B
	3/IEC		G321CL	G322CL	
	5/AHA		G511CL	G512CL	CN158A
	5/IEC		G521CL	G522CL	

ECG Leadwires/心电分线

Type Designation 型号类型	Leads/label 导联数	Trunk cable Conn. 分线盒接头	Pinch 夹式	Snap 扣式	Leadwires Conn. 导联线接头
GE-Datex Ohmeda 德恩	3/AHA	CN163	G311DX	G312DX	CN165
	3/IEC		G321DX	G322DX	
	5/AHA	CN164	G511DX	G512DX	CN166
	5/IEC		G521DX	G522DX	
Philips-HP 飞利浦	3/AHA	CN173	G311AA	G312AA	CN175
	3/IEC		G321AA	G322AA	
	5/AHA	CN174	G511AA	G512AA	CN176
	5/IEC		G521AA	G522AA	
	3/AHA	CN179	G311PH	G312PH	CN181
	3/IEC		G321PH	G322PH	
	5/AHA	CN180	G511PH	G512PH	CN182
	5/IEC		G521PH	G522PH	

ECG Leadwires

Type Designation	Leads/label	Trunk cable Conn.	Pinch	Snap	Leadwires Conn.
Philips-HP	10/AHA	 CN600	DMS511PV	DMS512PV	 CN188
	10/IEC		DMS521PV	DMS522PV	
Mindray-Datascope	3/AHA	 CN184	DMS311DT	DMS312DT	 CN168-3
	3/IEC		DMS321DT	DMS322DT	
	5/AHA		DMS511DT	DMS512DT	
	5/IEC		DMS521DT	DMS522DT	
	3/AHA	 CN159	DMS311MD	DMS312MD	 CN161
	3/IEC		DMS321MD	DMS322MD	 CN161
	5/AHA	 CN160	DMS511MD	DMS512MD	 CN162
	5/IEC		DMS521MD	DMS522MD	 CN162
5/AHA	 CN160-SR	DMS511MD -FM	DMS512MD -FM	 CN162	
5/IEC		DMS521MD -FM	DMS522MD -FM		

ECG Leadwires

Type Designation	Leads/label	Trunk cable Conn.	Pinch	Snap	Leadwires Conn.
GE-Marquette	10/AHA	 CN251	DMS111MQ	DMS121MQ	 CN168
	10/IEC		DMS112MQ	DMS122MQ	 CN168
	3/AHA	 CN167	DMS311MQ	DMS312MQ	 CN168-3
	3/IEC		DMS321MQ	DMS322MQ	 CN168-3
	5/AHA	 CN167	DMS511MQ	DMS512MQ	 CN168
	5/IEC		DMS521MQ	DMS522MQ	 CN168
	6/AHA	 CN136	DMS611MQ	DMS612MQ	 CN137
	6/IEC		DMS621MQ	DMS622MQ	 CN137
	10/AHA	 CN251	DMS111MQ	DMS112MQ	 CN168
	10/IEC		DMS121MQ	DMS122MQ	 CN168
	10/AHA	 CN221	DMS111GE	DMS112GE	 CN168B
	10/IEC		DMS121GE	DMS122GE	 CN168B

ECG



ECG Leadwires

Type Designation	Leads/label	Trunk cable Conn.	Pinch	Snap	Leadwires Conn.
<b>Nihon Kohden</b>	3/AHA	 CN171	DMS311NK	DMS312NK	 CN172-3
	3/IEC		DMS321NK	DMS322NK	
	5/AHA		DMS611NK	DMS612NK	 CN172
	5/IEC		DMS621NK	DMS622NK	
	3/AHA	 CN169	DMS311NH	DMS312NH	 CN170-3
	3/IEC		DMS321NH	DMS322NH	
	5/AHA		DMS511NH	DMS512NH	 CN170
	5/IEC		DMS521NH	DMS522NH	
<b>Spacelabs</b>	3/AHA	 CN177	DMS311SL	DMS312SL	 CN178
	3/IEC		DMS321SL	DMS322SL	 CN178
	5/AHA		DMS511SL	DMS512SL	 CN178
	5/IEC		DMS521SL	DMS522SL	 CN178































ECG Leadwires

Type Designation	Leads/label	Trunk cable Conn.	Pinch	Snap	Leadwires Conn.
<b>Drager-Siemens</b>	3/AHA	 CN159	DMS311DR	DMS312DR	 CN161
	3/IEC		DMS321DR	DMS322DR	 CN161
	5/AHA	 CN160	DMS611DR 5956490	DMS612DR	 CN162
	5/IEC		DMS621DR 5956508	DMS622DR	 CN162

ECG

	3/AHA	5/AHA	10/AHA	3/IEC	5/IEC	10/IEC
<b>Pinch</b>						
<b>Snap</b>						

Neonate ECG Leadwires

Type Designation	Leads/label	Trunk cable Conn.	Mini Pinch	Leadwires Conn.	Photo
Philips	3/AHA		DMS315AA	 CN175	
	3/IEC	CN173	DMS325AA	 CN175	
	5/AHA		DMS515AA	 CN176	
	5/IEC	CN174	DMS525AA	 CN176	
Din type	3/AHA		DMS315DN	 CN153	
	3/IEC	CN151	DMS325DN	 CN153	
	5/AHA		DMS515DN	 CN153	
	5/IEC	CN152	DMS525DN	 CN153	
Drager-Siemens	3/AHA		DMS315DR	 CN161	
	3/IEC	CN159	DMS325DR	 CN161	
	5/AHA		DMS515DR	 CN162	
	5/IEC	CN160	DMS525DR	 CN162	

Neonate ECG Leadwires


Type Designation	Leads/label	Trunk cable Conn.	Mini Pinch	Leadwires Conn.	Photo
Mindray	3/AHA		DMS315MD	 CN161	
	3/IEC	CN159	DMS325MD	 CN161	
	5/AHA		DMS515MD	 CN162	
	5/IEC	CN160	DMS525MD	 CN162	

ECG

Multi-Parameter Patient Cable

P/N	Description	Used with the OEM and models	Pictures
DMS6136TR 3368391 5950196	Drager-Siemens Multimed 5Pod patient cable ,2.5M5-leadECG, VS/Siemens Spo2/Temp/Resp AHA&IEC	Siemens/Drager Infinity6000/ 7000/9000 Drager Infinity Gamma/Gamma XL/Gamma X XL Delta/DeltaXL Kappa/ KappaXL T Vista/VistaXL	
DMS6136SR 5590539	Drager-Siemens Neomed patient cable 2.5M ECG/Resp/Spo2 Dual Temp/FiO2 AHA&IEC	Drager Infinity Gamma/ GammaXL/GammaX XL Delt/DelteXL Kappa/ KappaXL T Vista/VistaXL	
DMS6136DR MS20093	Drager-Siemens Multimed Plus patient cable 2.5M Spo2/Temp/6-lead ECG/Resp AHA&IEC	Drager Infinity Gamma/ GammaXL/Gamma X XL Delta/DeltaXL Kappa/ KappaXL T Vista/VistaXL	
DMS321DR-6P MS14556	Drager-Siemens single pin multi-lead wire,3leads Pinch (Snap available) IEC(AHA available)	Multimed Plus MS20093	
DMS521DR-6P MS14560	Drager-Siemens single pin multi-lead wire,5leads Pinch(Snap available) IEC(AHA available)		
DMS621DR MS14682	Drager-Siemens single pin multi-lead wire,6leads Pinch(Snap available) IEC(AHA available)		

Holter Recorder ECG

Models	Connector	Leads/label	P/N	OEM #
Normal DIN 1.5	 CN153	7/AHA Leadwire-Snap	DMS712DN	
		10/AHA Leadwire-Snap	DMS112DN	
		10/IEC Leadwire-Snap	DMS122DN	
Contec	 CN153	5/AHA Leadwire-Snap	DMS5193S	
		10/IEC Fixed-Snap	DMS12123S	
Mortala	 CN122	10/AHA Fixed-Snap	DMS1180S	
		10/IEC Fixed-Snap	DMS1280S	
Mortala H3		5/AHA Fixed-Snap	DMS5190S	
		5/IEC Fixed-Snap	DMS5290S	
Mindray	 CN434	5/AHA Fixed-Snap	DMS5183S	
Landcom	 CM10	10/AHA Fixed-Snap	DMS1182S	
		10/IEC Fixed-Snap	DMS1282S	
Philips Digitrak XT	 CN134	5/AHA Fixed-Snap	DMS5188S	989803157491 (Medium)
		5/IEC Fixed-Snap	DMS5288S	
Philips Digitrak Plus	 CN255	5/AHA Fixed-Snap	DMS5189S	M4725A
		5/IEC Fixed-Snap	DMS5289S	
Philips ST80I		10/AHA Fixed-Snap	DMS11148S	
Philips		5/IEC Grabber ECG+SpO2	DMS0040	989803171951
		6/IEC Grabber ECG+SpO2	DMS0040B	

ECG

Holter Recorder ECG

Models	Connector	Leads/label	P/N	OEM #
Philips		5/AHA Grabber ECG+SpO2	DMS0040C	989803171851
GE Seer Light	 CN256	5/AHA Fixed-Snap	DMS5191S	2008594-001
		7/AHA Fixed-Snap	DMS7191S	2008594-002
		7/IEC Fixed-Snap	DMS7291S	2008594-004
		10/AHA Fixed-Snap	DMS1191S	
JincoMed	 CP14	10/AHA Fixed-Snap	DMS1181S	
		10/IEC Fixed-Snap	DMS1281S	
BI(Biomedical Instruments) BI-9000TL	 CP14	7/AHA Fixed-Snap	DMS7187S	
		10/AHA Fixed-Snap	DMS1184S	
BI(Biomedical Instruments) BI9800/Edan Advanced HT-1000	 CN428	5/AHA Fixed-Snap	DMS5195S	
		7/AHA Fixed-Snap	DMS7195S	
		10/AHA Fixed-Snap	DMS1195S	
DMS 300-3A	 CN427	5/AHA Fixed-Snap	DMS5185S	
		7/AHA Fixed-Snap	DMS7185S	
		10/AHA Fixed-Snap	DMS1185S	
DMS 300-3A DMS 300-2	 CN430	5/AHA Fixed-Snap	DMS51113S	
		7/AHA Fixed-Snap	DMS71113S	
		10/AHA Fixed-Snap	DMS1113S	
Schiller MT-101, MT-200	 CN429	4AHA/IEC Fixed-Snap	DMS422SC	
		6AHA/IEC Fixed-Snap	DMS622SC	
Schiller MS-12	 CN472	10AHA/IEC Fixed-Snap	DMS11147S	
Borsam BS6930	 CN465	5/AHA Fixed-Snap	DMS1129S	

Holter Recorder ECG

Models	Connector	Leads/label	P/N	OEM #
Borsam iTengo+	 CN466	5/AHA Fixed-Snap	DMS51130S	
Biocare IH-12PLUS	 CN467	10/AHA Fixed-Snap	DMS11137S	
		10/IEC Fixed-Snap	DMS12137S	
BTL uk BTL-08	 CN468	10/AHA Fixed-Snap	DMS51138S	
		10/IEC Fixed-Snap	DMS52138S	
Carewell	 CN469	10/AHA Fixed-Snap	DMS11142S	
		10/IEC Fixed-Snap	DMS12142S	
Edan SE-2003 ,SE-2012	 CN470	5/AHA Fixed-Snap	DMS1125S	
		7/AHA Fixed-Snap	DMS71125S	
Edan SE-1010	 CN471	10/IEC Fixed-Snap	DMS12140S	
Shanghai Quntian PI200D-BY	 CN459+CN460	10/IEC Fixed-Snap	DMS12146S	
Fukuda Denshi LX-2120	 CN473	4/AHA Fixed-Snap	DMS412FT	
		4/IEC Fixed-Snap	DMS422FT	

ECG

EEG Cables

MEDKE P/N	Description	Package	Photo
DMS0001	Reusable grey EEG cup electrodes, 1.5M Golden plated Cup electrode EEG cord	1 pack/set of 10lead	
DMS0001-B	Reusable color-coded EEG electrodes, 1.5M Golden plated Cup electrode EEG cord	1 pack/set of 10lead	






## EKG Cables & Accessories

## EKG Cables & Accessories

- Fixed One-piece, Trunk cable, Leadwire types available
- Latex-free
- Banana, Needle, Snap, Pinch available
- AHA and IEC color coding available
- Different types are Compatible with most of patient monitors
- Easy to clean and repair



ECG

Product types	Photo	P/N Example	P/N Guide
Fixed one-piece EKG cables		DMS1113B	K: EKG cables 1: 10 leads EKG cable 1: AHA (2: IEC) 13: Philips (01: Burdick) B: Banana (N: Needle)
Trunk cable		DMS1106MQ	K: EKG cables 1: 10 leads EKG cable 1: AHA (2: IEC) 06: GE-Marquette MQ: MQ type
Leadwire		DMS114MQ	K: EKG cables 1: 10 leads EKG cable 1: AHA (2:IEC) 4: Banana (3:Needle,2:Snap,1:Pinch) MQ: MQ type
Warranty: one-year warranty Package: 1pcs/bag			

## Disposable IBP Transducer

- Single patient use
- Easy to use with interface cables
- Abbott/Medex, Edward, Philips/ B.Bruan, B.D, Utah types available
- Latex free











Photos	P/N	Description	Package	MOQ
	DMS-101	Abbott/Medex pressure transducer kits	1pcs/pack	10pcs
	BB-102	Edward pressure transducer kits	1pcs/pack	10pcs
	DMS-103	Philips/B.Bruan pressure transducer kits	1pcs/pack	10pcs
	DMS-104	B.D pressure transducer kits	1pcs/pack	10pcs
	DMS-105	Utah pressure transducer kits	1pcs/pack	10pcs



## Temperature Probes










EKG Accessories

	10/AHA	10/IEC
Pinch		
Snap		
Needle		
Banana		

EKG Accessories

P/N	Description	Color	Package	Pictures 片
DMS9201	Adult limb clamps electrode with screw, Ag/AgCl, universal ,3mm/4mm, IEC GE OEM: 2029894-001	Red/black/ green/yellow	4 pcs/pack	

P/N	Description	Color	Package	Pictures
DMS9203	Adult limb clamps electrode, Nickel, 3mm/4mm, IEC	Red/black/ green/yellow	4 pcs/pack	
DMS9205	Adult limb clamp electrode with Nickel Universal , metal spring	Red/black/ green/yellow	4 pcs/pack	
DMS7201	Pediatric limb clamps electrode with screw, Ag/AgCl, universal ,3mm/4mm, IEC	Red/black/ green/yellow	4 pcs/pack	
DMS9202	Adult suction electrode with screw, Ag/AgCl, universal, 3mm/4mm Bulb diameter: 28mm GE OEM: 2029892-001	Blue	6 pcs/pack	
DMS9204	Adult suction electrode, Nickel, 3mm/4mm Bulb diameter: 28mm GE OEM: 21714401	Blue	6 pcs/pack	
DMS9207	Adult suction electrode , Nickel universal red,yellow,purple,green,black, brown Bulb diameter: 28mm	Red/yellow/purple/ green/black/brown	6 pcs/pack	
DMS7207	Pediatric suction electrode with screw, Ag/AgCl, universal, 3mm/4mm Bulb diameter: 21mm	Blue	6 pcs/pack	
DMS9208	Adult suction electrode, Nickel, Bulb diameter: 28mm	Blue	6 pcs/pack	

ECG

Disposable ECG

P/N	Description	Package	MOQ	Photo
DMS-101	Disposable adult foam electrode, Solid gel Diameter: 55mm	50pcs/bag	1 pack	
DMS-102	Disposable pediatric foam electrode, Solid gel Diameter: 43mm	50pcs/bag	1 bag	
DMS-103	Disposable infant foam electrode, Solid Diameter: 30mm	100pcs/bag	1 bag	
DMS0004	Disposable ECG patch DIN 1.5mm 60cm in length , Dia: 22mm Red/White/black ,AHA	3pcs/bag, 1set / bag	200 set	
DMS0044B	Disposable ECG patch DIN 1.5mm 60cm in length , Dia: 22mm Red/yellow/green ,IEC	3pcs/bag, 1set / bag	200 set	
DMS0004C	Disposable ECG patch DIN 1.5mm 60cm in length , Dia: 22mm Red*2/yellow*2/green ,IEC	5pcs/bag, 1set / bag	200 set	

ECG Adapter

P/N	Description	Package	Pictures
DMS01	Vet ECG adapter use with snap/grabber end	5pcs/bag	
DMS02	Universal vet ECG Adapter Use with snap/Grabber/3mm needle/4mm banana	5pcs/bag	
DMS05	ECG adapter clip 4mm Banana to Grabber/Pinch end	10pcs/bag	
DMS05-B	ECG adapter clip 3mm needle to Grabber/Pinch end	10pcs/bag	
DMS12	ECG adapter clip 4mm Banana to Snap end	10pcs/bag	
DMS12-C	ECG Adapter Clips 3mm needle	10pcs/bag	
DMS15	Multifunction Electrodes Adapter Locks on to Tab electrodes Fits Banana plugs (4mm socket)	10pcs/bag	
DMS16	Multifunction Electrodes Adapter Locks on to Tab electrodes Fits Needle plugs (3mm socket)	10pcs/bag	

ECG

ECG Adapter

P/N	Description	Package	Pictures
DMS04	Alligator Clip /red adapter 4mm Banana	10pcs/bag	
DMS13	ECG Adapter Clips 4mm Banana	10pcs/bag	
DMS18	Alligator Clip /blue adapter 3mm needle	10pcs/bag	
DMS30NS	ECG Adapter wire 3mm Needle to Snap	10pcs/bag	
DMS40BS	ECG/EKG Adapter Cable Banana to Snap for Disposable ECG Cable	10pcs/bag	
DMS30NP	ECG Adapter wire 3mm Needle to Grabber/Pinch	10pcs/bag	
DMS40BP	ECG/EKG Adapter Cable Banana to Grabber/pinch for Disposable ECG Cable	10pcs/bag	



Reusable NIBP Cuffs

## NIBP Cuffs, Connectors And Adapter Hoses

- Reusable and disposable cuffs are available with Single and double tube
- Easy to clean
- Latex free



Product types	Photo	P/N Example	P/N Guide
NIBP disposable cuffs		DMS0105	C: NIBP cuffs 0: Disposable cuff 1: single tube (2: two tubes) 05: cuff size
NIBP reusable cuffs with bag		DMS1711	C: NIBP cuffs 1: Reusable cuff nylon, 2 Reusable cuff PU brown, 6 Reusable cuff PU blue 7: Adult (6: Pediatric; 5: Infant; 3: Neonate; 8: LargeAdult; 9: Thigh) 1: single tube (2: two tubes) 1: with bag
NIBP reusable cuffs without bag		DMS1710	C: NIBP cuffs 1: Reusable cuff brown 7: Adult (6: Pediatric; 5: Infant; 3: Neonate; 8: LargeAdult; 9: Thigh) 1: single tube (2: two tubes) 0: without bag
NIBP adapter hoses		DMS3907D	H: NIBP adapter hoses 39: CB39 connector 07: CB07 connector D: Double tubes
Warranty: six-month warranty for NIBP cuffs; six-month warranty for NIBP Adapter Hoses; Package: 1pcs/bag			

## Disposable NIBP Cuffs

Size	Limb Circumference	Single tube P/N	Double tube P/N	MOQ	Image
Neonate 1	3-6cm	DMS0101-44	DMS0201-4444	20pcs	
Neonate 2	4-8cm	DMS0102-44	DMS0202-4444	20pcs	
Neonate 3	6-11cm	DMS0103-44	DMS0203-4444	20pcs	
Neonate 4	7-13cm	DMS0104-44	DMS0204-4444	20pcs	
Neonate 5	8-15cm	DMS0105-44	DMS0205-4444	20pcs	
neonate 1	3-6cm	DMS0301-50	DMS0401-5050	20pcs	
neonate 2	4-8cm	DMS0302-50	DMS0402-5050	20pcs	
neonate 3	6-11cm	DMS0303-50	DMS0403-5050	20pcs	
neonate 4	7-13cm	DMS0304-50	DMS0404-5050	20pcs	
neonate 5	8-15cm	DMS0303-44	DMS0403-44	20pcs	
Infant	9.8-13.3cm	C0107-43/C0107-34/C0107-32	C0207-4343/C0207-3434/C0207-3232	5pcs	
Small Child	12.4-16.8cm	C0108-43/C0108-34/C0108-32	C0208-4343/C0208-3434/C0208-3232	5pcs	
Child	15.8-21.3cm	C0109-43/C0109-34/C0109-32	C0209-4343/C0209-3434/C0209-3232	5pcs	
Small adult	20-27cm	C0110-43/C0110-34/C0110-32	C0210-4343/C0210-3434/C0210-3232	5pcs	
Adult	25.3-34.3cm	C0111-43/C0111-34/C0111-32	C0211-4343/C0211-3434/C0211-3232	5pcs	
Large Adult	32.1-43.4cm	C0112-43/C0112-34/C0112-32	C0212-4343/C0212-3434/C0212-3232	5pcs	
Thigh	46-66cm	C0113-43/C0113-34/C0113-32	C0213-4343/C0213-3434/C0213-3232	5pcs	
Please match right plastic NIBP Connector for cuffs					

## Reusable D-ring Cuff For Digital Sphygmomanometer ABPM

Size	Limb Circumference	Single tube P/N	MOQ	Image
Pediatric	18-26cm	DMS8611F	1pcs/bag	
Adult	22-32cm	DMS8711F	1pcs/bag	
Large Adult	32-43cm	DMS8811F	1pcs/bag	
Large Adult Long	36-45cm	DMS8911F	1pcs/bag	

NIBP Cuffs

P/N	Description (Bladder)	P/N	Description (Bladderless)	With connector
DMS1721-3636	Adult, Nylon, single tube, 25-35cm	DMS1720-3636	Adult, TPU&Nylon, double tube, 25.3-34.3cm	 
DMS1621-3636	Pediatric, Nylon, double tube, 18-26cm	DMS1620-3636	Small adult, TPU&Nylon, double tube, 20-27cm	
DMS1521-3636	Infant, Nylon, double tube, 10-19cm	DMS1520-3636	Infant, TPU&Nylon, double tube, 12.4-16.8cm	
		DMS1420-3636	Child, TPU&Nylon, double tube, 15.8-21.3cm	
DMS1321-3636	Neonate, Nylon, double tube, 6-11cm	DMS1320-3636	Neonate, TPU&Nylon, double tube, 6.9-11.7cm	
DMS1621-3636	Large Adult, Nylon, double tube, 33-47cm	DMS1620-3636	Large Adult, TPU&Nylon, double tube, 32.1-43.4cm	
DMS1921-3636	Thigh, Nylon, double tube, 46-66cm	DMS1920-3636	Thigh, TPU&Nylon, double tube, 46-66cm	
P/N	Description (Bladder)	P/N	Description (Bladderless)	With connector
DMS1721-3232	Adult, Nylon, double tube, 25-35cm	DMS1720-3232	Adult, TPU&Nylon, double tube, 25.3-34.3cm	 
DMS1621-3232	Pediatric, Nylon, double tube, 18-26cm	DMS1620-3232	Small adult, TPU&Nylon, double tube, 20-27cm	
DMS1721-3232	Infant, Nylon, double tube, 10-19cm	DMS1720-3232	Infant, TPU&Nylon, double tube, 12.4-16.8cm	
		DMS1420-3232	Child, TPU&Nylon, double tube, 15.8-21.3cm	
DMS1321-3232	Neonate, Nylon, double tube, 6-11cm	DMS1320-3232	Neonate, TPU&Nylon, double tube, 6.9-11.7cm	
DMS1821-3232	Large Adult, Nylon, double tube, 33-47cm	DMS1820-3232	Large Adult, TPU&Nylon, double tube, 32.1-43.4cm	
DMS1921-3232	Thigh, Nylon, double tube, 46-66cm	DMS1920-3232	Thigh, TPU&Nylon, double tube, 46-66cm	
P/N	Description (Bladder)	P/N	Description (Bladderless)	With connector
DMS1721-3132	Adult, Nylon, double tube, 25-35cm	DMS1720-3132	Adult, TPU&Nylon, double tube, 25.3-34.3cm	 
DMS1621-3132	Pediatric, Nylon, double tube, 18-26cm	DMS1620-3132	Small adult, TPU&Nylon, double tube, 20-27cm	
DMS1521-3132	Infant, Nylon, double tube, 10-19cm	DMS1520-3132	Infant, TPU&Nylon, double tube, 12.4-16.8cm	
		C1420-3132	Child, TPU&Nylon, double tube, 15.8-21.3cm	
DMS1321-3132	Neonate, Nylon, double tube, 6-11cm	DMS1320-3132	Neonate, TPU&Nylon, double tube, 6.9-11.7cm	
DMS1821-3132	Large Adult, Nylon, double tube, 33-47cm	DMS1820-3132	Large Adult, TPU&Nylon, double tube, 32.1-43.4cm	
DMS1921-3132	Thigh, Nylon, double tube, 46-66cm	DMS1920-3132	Thigh, TPU&Nylon, double tube, 46-66cm	

NIBP Cuffs

P/N	Description (Bladder)	P/N	Description (Bladderless)	With connector
DMS1711-43	Adult, Nylon, single tube, 25-35cm	DMS1710-43	Adult, TPU&Nylon, single tube, 25.3-34.3cm	
DMS1611-43	Pediatric, Nylon, single tube, 18-26cm	DMS1610-43	Small adult, TPU&Nylon, single tube, 20-27cm	
DMS1511-43	Infant, Nylon, single tube, 10-19cm	DMS1510-43	Infant, TPU&Nylon, single tube, 12.4-16.8cm	
		DMS1410-43	Child, TPU&Nylon, single tube, 15.8-21.3cm	
DMS1311-43	Neonate, Nylon, single tube, 6-11cm	DMS1310-43	Neonate, TPU&Nylon, single tube, 6.9-11.7cm	
DMS1811-43	Large Adult, Nylon, single tube, 33-47cm	DMS1810-43	Large Adult, TPU&Nylon, single tube, 32.1-43.4cm	
DMS1911-43	Thigh, Nylon, single tube, 46-66cm	DMS1910-43	Thigh, TPU&Nylon, single tube, 46-66cm	
P/N	Description (Bladder)	P/N	Description (Bladderless)	With connector
DMS1711-33	Adult, Nylon, single tube, 25-35cm	DMS1710-33	Adult, TPU&Nylon, single tube, 25.3-34.3cm	
DMS1611-33	Pediatric, Nylon, single tube, 18-26cm	DMS1610-33	Small adult, TPU&Nylon, single tube, 20-27cm	
DMS1511-33	Infant, Nylon, single tube, 10-19cm	DMS1510-33	Infant, TPU&Nylon, single tube, 12.4-16.8cm	
		DMS1410-33	Child, TPU&Nylon, single tube, 15.8-21.3cm	
DMS1311-33	Neonate, Nylon, single tube, 6-11cm	DMS1310-33	Neonate, TPU&Nylon, single tube, 6.9-11.7cm	
DMS1811-33	Large Adult, Nylon, single tube, 33-47cm	DMS1810-33	Large Adult, TPU&Nylon, single tube, 32.1-43.4cm	
DMS1911-33	Thigh, Nylon, single tube, 46-66cm	DMS1910-33	Thigh, TPU&Nylon, single tube, 46-66cm	
P/N	Description (Bladder)	P/N	Description (Bladderless)	With connector
DMS1711-32	Adult, Nylon, single tube, 25-35cm	DMS1710-32	Adult, TPU&Nylon, single tube, 25.3-34.3cm	
DMS1611-32	Pediatric, Nylon, single tube, 18-26cm	DMS1610-32	Small adult, TPU&Nylon, single tube, 20-27cm	
DMS1511-32	Infant, Nylon, single tube, 10-19cm	DMS1510-32	Infant, TPU&Nylon, single tube, 12.4-16.8cm	
		C1410-32	Child, TPU&Nylon, double tube, 15.8-21.3cm	
DMS1311-32	Neonate, Nylon, single tube, 6-11cm	DMS1310-32	Neonate, TPU&Nylon, single tube, 6.9-11.7cm	
DMS1811-32	Large Adult, Nylon, single tube, 33-47cm	DMS1810-32	Large Adult, TPU&Nylon, single tube, 32.1-43.4cm	
DMS1911-32	Thigh, Nylon, single tube, 46-66cm	DMS1910-32	Thigh, TPU&Nylon, single tube, 46-66cm	

NIBP

NIBP Connectors

P/N	DMS01	DMS01	DMS03	DMS04	DMS05	DMS06	DMS07
Photo							
Match	CB01->CB02/CN44	CB02->CB01	CB03->CB04	CB04->CB03	CB05->CB06	CB06->CB05	CB07->CB36/CB16
Use with	For Neonate air hose conn.	CAS, Drager, Kontron, Physio control, Schiller /Bruker/Odam	Mennen (Old), Schiller/Bruker/Odam, Spacelabs, Nihon Kohden Monitor plug	BCI, Datascope(Old)	Monitor plug	Colin, Gambro, Tyco( NPB 5500)	Critikon, Invivo, Kontron, Nihon Kohden GE:300619
P/N	DMS08	DMS09	DMS10	DMS11	DMS12	DMS13	DMS14
Photo							
Match			CB10->CB11	CB11->CB10			CB14->CB70
Use with	ID Barbed to ID Barbed, Metal	Philips/HP Monitor plug	GE:330060 Mindray, Schiller	Philips, Siemens, Spacelabs, CSI, Welch Allyn, Datascope GE:330059 WA:5082-179	Drager Monitor plug	Monitor plug	PACE TECH Monitor plug
P/N	DMS16	DMS30	DMS31	DMS32	DMS33	DMS34	DMS36
Photo							
Match	CB16->CB07		CB31->CB32	CB32->CB31	CB33->CB34	CB33->CB34	CB36->CB07
Use with	Critikon, Invivo, Kontron, Mennen, Nihon Kohden, Welch Allyn GE:300664 WA:5082-164	ID Barbed to ID Barbed, Plastic GE:300670	Welch Allyn, Marquette Datex-Ohmeda, Nihon Kohden GE:330091 WA:5082-178	Welch Allyn Marquette, Datex, Zoll GE:330064 WA:5082-172	BCI, Datascope GE:300665 WA:5082-165	BCI, Datascope GE:300668 WA:5082-168	Critikon, Invivo, Kontron, Mennen, Nihon Kohden, Welch Allyn GE:300664 WA:5082-164
P/N	DMS37	DMS38	DMS39	DMS40	DMS41	DMS68	DMS42
Photo							
Match							
Use with	GE AS/3, CS/3, S/5, Eagle4000, Cardiocap 5 Monitor plug	GE-Dinamap Monitor plug	GE-Dash, Eagle, Solar Monitor plug	Nihon Kohden For adult/ pediatric Monitor plug	Nihon Kohden For neonate/ infant Monitor plug	Nihon Kohden visom/SVM-7523	for neonate double tube disposable BP cuffs

NIBP Connectors

P/N	DMS43	DMS44	DMS45	DMS46	DMS49	DMS50	DMS51
Photo							
Match	CB43->CB10	CB44->CB01		CB46->CB70		CB50->CB01	CB51->CB52
Use with	Philips, Siemens, Spacelabs, CSI, Datascope WA:5082-184	Neonate disposable cuff WA:5082-176	Use with all NIBP connector	Pace Tech	Universal Omron connector	Neonate disposable cuff	GE Dinacllick
P/N	DMS57	DMS52	DMS55	DMS56	DMS59	DMS60	DMS62
Photo							
Match	CB57->CB52	CB52->CB51/ CB57		CB56-> CB66	CB59-> CB60	CB59-> CB60	
Use with	GE Dinacllick	GE new Dinacllick NIBP	Welch Allyn Monitor plug	Drager (Cuff side)	For Philips disposable neonate tube side	For Philips disposable neonate cuff side	For Welch Allyn cuff side
P/N	DMS63	DMS66	DMS70	DMS10-A	DMS11	DMS08-A	
Photo							
Match		CB56-> CB66	CB46 >CB70				
Use with	For Welch Allyn cuff side	Drager (Tube side)	M&B	(OEM)SAIRAN	(OEM)SAIRAN	(OEM)SAIRAN	

NIBP



## NIBP Adapter Hoses

### NIBP Adapter Hoses

Models	Tube	Monitor side Conn.	Cuff side Conn.	Options	P/N	OEM #
Air shields Athena	2	 CB11	 CB07	P/N for Adult&Pedia.	DMS1107D	
			 CB42	P/N for Infant&Neo.	DMS142ND	
Colin 8800	1	 CB05	 CB06	P/N for Adult&Pedia.	DMS0506S	
			 CB01	P/N for Infant&Neo.	DMS0501NS	
CSI 507,1100	2	 CB37	 CB07	P/N for Adult&Pedia.	DMS3707D	558
			 CB42	P/N for Infant&Neo.	DMS3742ND	
CSI	1	 CB05	 CB10	P/N for Adult&Pedia.	DMS0510S	709
			 CB01	P/N for Infant&Neo.	DMS0501NS	
CSI	1	 CB10	 CB10	P/N for Adult&Pedia.	DMS1010S	706
			 CB01	P/N for Infant&Neo.	DMS1001NS	
CSI 507E,508, Scholar II (See Colin)	1	 CB05	 CB06	P/N for Adult&Pedia.	DMS0506S	
			 CB01	P/N for Infant&Neo.	DMS0501NS	

NIBP

NIBP Adapter Hoses

Models	Tube	Monitor side Conn.	Cuff side Conn.	Options	P/N	OEM #
Drager-Siemens SC6000-7000/ 8000/9000 Infinity	1		CB10 	P/N for Adult&Pedia.	DMS1210S	1275275 EH50U
			CB01 	P/N for Infant&Neo.	DMS1201NS	2870298
GE-Critikon Dinamap 1846SX, Dinamap8100, Dinamap 9300XL	2		CB07 	P/N for Adult&Pedia.	DMS0707D	008841 E9050JY
			CB42 	P/N for Infant&Neo.	DMS0742ND	008840
GE-Critikon (See CSI) GE-Datex Ohmeda  AS/3,CS/3, Cardiicap 5, Cardiicap II, Light Monitor S/5	2		CB07 	P/N for Adult&Pedia.	DMS3707D	9461-214
			CB42 	P/N for Infant&Neo.	DMS3742ND	877514
GE-Datex Ohmeda  AS/3,CS/3, Cardiicap 5, Cardiicap II, Light Monitor S/5	2		CB31  CB32 	P/N for Adult&Pedia.	DMS37312D	877235558
GE-Marquette Eagle 4000, Tram 100/200/ 250/300/350, Corometrics 120	2		CB07 	P/N for Adult&Pedia.	DMS3907D	9461-203 E9001TE
			CB42 	P/N for Infant&Neo.	DMS3942ND	877514
GE-Marquette Eagle 3000, DASH 2000/3000/4000, Tram 400/450/451, Tram 800/850/851, Corometrics 120	2		CB31 	P/N for Adult&Pedia.	DMS391D	414873-001 E9001TC 2017008-001
			CB42 	P/N for Infant&Neo.	DMS3942ND	414874-001 2017009-001

NIBP Adapter Hoses

Models	Tube	Monitor side Conn.	Cuff side Conn.	Options	P/N	OEM #
GE-Marquette Dinamap MPS Series,Dinamap PRO100/110/200/ 210/300/310/400 Dinamap ProCare 100/120/200/220/ 300/320/400	2		CB07 	P/N for Adult&Pedia.	DMS3807D	107365
			CB42 	P/N for Infant&Neo.	DMS3842ND	107368
			CB51 	P/N for Adult&Pedia.	DMS3851D	2058205-001
Invivo, MDE 1400,5600(See Nihon kohden)	2		CB07 	P/N for Adult&Pedia.	DMS0707D	
			CB42 	P/N for Infant&Neo.	DMS0742ND	
Kontron (AVCO)  Minimon 7137B, Minimon 7137Plus (See Critikon )	2		CB07 	P/N for Adult&Pedia.	DMS0707D	
			CB42 	P/N for Infant&Neo.	DMS0742ND	
Medtronic-Physio Contro Lifepak 12	1		CB10 	P/N for Adult&Pedia.	DMS1110S	
			CB01 	P/N for Infant&Neo.	DMS1101NS	
Mindray PM-600,PM-6000, PM-7000,PM-8000, PM-8000 Express, PM-9000 Express, VS-800(See CSI) Mindray-Datascope AccuTorr Plus,Duo, Passport 2,Passport 2LT,Spectrum, Spectrum OR Trio	1		CB10 	P/N for Adult&Pedia.	DMS1010S	M3918A 0683-04-0003 0683-04-0004
			CB01 	P/N for Infant&Neo.	DMS1001NS	

NIBP



NIBP Adapter Hoses

Models	Tube	Monitor side Conn.	Cuff side Conn.	Options	P/N	OEM #
Mindray BeneView T5, BeneView T6, BeneView T8, MPM Module (See Philips)	1		CB10 	P/N for Adult&Pedia.	DMS0910S	
			CB01 	P/N for Infant&Neo.	DMS0901NS	
Mindray-Datascope 3000, Accutorr 1, Accutorr 2, Accutorr 3, Accutorr 4, AccutorrPlus, ExpertPassport, Passport 2, Passport2LT, Spectrum(See CSI )	1		CB10 	P/N for Adult&Pedia.	DMS0510S	0683-04-0006
			CB01 	P/N for Infant&Neo.	DMS0501NS	
Mindray-Datascope 2100,2200I, Accutorr 1, Accutorr2, Accutorr 3, Accutorr 4, AccutorrPlus, Passort 3L, Passport,Passport2, Passport 2LT, Passport 5L, Spectrum	1		CB10 	P/N for Adult&Pedia.	DMS0610S	0683-04-0007
			CB01 	P/N for Infant&Neo.	DMS0601NS	
Nellcor 3100,4000	1		CB10 	P/N for Adult&Pedia.	DMS1010S	
			CB01 	P/N for Infant&Neo.	DMS1001NS	
Nihon Kohden AP-860PA Module, BSM-2101,BSM- 2102,BSM-8301, BSM-8302,BSM- 8500,BSM-8502, BSM-8800, Lifescope 14, Lifescope 9 Lifescope L	2		CB07 	P/N for Adult&Pedia.	DMS0707D	6143-001321
			CB42 	P/N for Infant&Neo.	DMS0742ND	
Nihon Kohden AP-851PA Module, BSM-8300, BSM-8301, BSM-8302, Lifescope 9	1		CB04 	P/N for Adult&Pedia.	DMS0304S	110088
			CB01 	P/N for Infant&Neo.	DMS0301NS	

NIBP Adapter Hoses

Models	Tube	Monitor side Conn.	Cuff side Conn.	Options	P/N	OEM #
Nihon Kohden Cb40 is different from Cb41	2		CB32 	P/N for Adult&Pedia.	DMS4032D	YN-900P YN-901P
			CB01 	P/N for Infant&Neo.	DMS4101D	YN-920P YN-921P
Philips-HP 78352A,78352C, 78354A,78354C M1008A,M1008B, M3000	1		CB10 	P/N for Adult&Pedia.	DMS0910S	M1599B M1598B
			CB01 	P/N for Infant&Neo.	DMS0901NS	M1597B M1596B
Philips-HP Compact A1, Compact A3	1		CB10 	P/N for Adult&Pedia.	DMS1010S	M3918A
			CB01 	P/N for Infant&Neo.	DMS1001NS	
Spacelabs Medical Ultraview 90496, Ultraview 91496	1		CB10 	P/N for Adult&Pedia.	DMS1010S	714-0018-00
			CB01 	P/N for Infant&Neo.	DMS1001NS	
Spacelabs Medical Ultraview 90496, Ultraview 91496	1		CB34 	P/N for Adult&Pedia.	DMS1134S	714-0019-01
			CB01 	P/N for Infant&Neo.	DMS1101NS	
Spacelabs Medical All (dual tube) (See Nihon Kohden)	2		CB07 	P/N for Adult&Pedia.	DMS0707D	
			CB42 	P/N for Infant&Neo.	DMS0742ND	

NIBP

NIBP Adapter Hoses

Models	Tube	Monitor side Conn.	Cuff side Conn.	Options	P/N	OEM #
Spacelabs Medical All (single luer) except PCI (See Datascope)	1		 CB10	P/N for Adult&Pedia.	DMS0610S	
			 CB01	P/N for Infant&Neo.	DMS0601NS	
Welch Allyn Atlas, Protocol (See Nihon Kohden)	1		 CB04	P/N for Adult&Pedia.	DMS0304S	
			 CB01	P/N for Infant&Neo.	DMS0301NS	
Welch Allyn Spot Vital Signs IXi	2		 CB07	P/N for Adult&Pedia.	DMS5507D	10144
			 CB01	P/N for Infant&Neo.	DMS5501D	
	2		 CB31	P/N for Adult&Pedia.	DMS5531D	4500-32
			 CB63	P/N for Adult&Pedia.	DMS5563D	4500-34
Welch Allyn Encore 202EL, Encore 204EL, Encore 206EL, Propaq 102EL, Propaq 104EL, Propaq 106EL, Protocol	1		 CB31	P/N for Adult&Pedia.	DMS0731S	
			 CB01	P/N for Infant&Neo.	DMS0701NS	
Zoll Medical Masim Series (See Datascope)	1		 CB10	P/N for Adult&Pedia.	DMS0610S	
			 CB01	P/N for Infant&Neo.	DMS0601NS	

NIBP Adapter Hoses

Models	Tube	Monitor side Conn.	Cuff side Conn.	Options	P/N	OEM #
Comen, Creative, Zondan, General, Jeteem, HRD, GMI	1		 CB10	P/N for Adult&Pedia.	DMS1310S	
			 CB01	P/N for Infant&Neo.	DMS1310NS	
Fukuda Denshi	1		 CB10	P/N for Adult&Pedia.	DMS6910S	
			 CB01	P/N for Infant&Neo.	DMS6910NS	

NIBP



## IBP Adapter Cables

- Compatible with most of monitors
- Latex free

### B 03 01 :

**B:** IBP cable;

**03:** Edward(01:B.Bruan; 02:BD; 04:Medex/Abbott; 05:Utah; 06:PVB 07:Argon; 08:Medex Logical; 09:USB; 10: Mindray) 09:USB; 10: mindray)

**01:** 6Pin General(02:Mindray-Datascope; 03:Drager)

**Warranty:**six-month warranty for IBP adapter cables; **Package:** 1pcs/bag



IBP

## IBP Cables



B. Bruan



BD



Edward



Medex/Abbott



Utah



PVB



Argon



Medex Logical









USB







Mindray

IBP Cables

Monitor	Options										Conn.
	B.Bruan	BD	Edward	Medex/Abbott	Utah	PVB	Argon	Medex Logical	USB	Mindray	
<b>Spacelabs</b>											
	B0116	B0216	B0316	B0416	B0516	B0616	B0716	B0816	B0916	B1016	 CN312
<b>Fukuda Denshi</b>											
	B0117	B0217	B0317	B0417	B0517	B0617	B0717	B0817	B0917	B1017	 CN123B
<b>CSI/Criticare</b>											
	B0119	B0219	B0319	B0419	B0519	B0619	B0719	B0819	B0919	B1019	 CP5-50R
<b>MKE 6P</b>											
	B0121	B0221	B0321	B0421	B0521	B0621	B0721	B0821	B0921	B1021	 CN101
<b>Biolight</b>											
	B0123	B0223	B0323	B0423	B0523	B0623	B0723	B0823	B0923	B1023	 CN314
<b>Comen (12pin)</b>											
	B0127	B0227	B0327	B0427	B0527	B0627	B0727	B0827	B0927	B1027	 CN315
<b>Bistos</b>											
	B0128	B0228	B0328	B0428	B0528	B0628	B0728	B0828	B0928	B1028	 CN316
<b>Edan</b>											
	B0129	B0229	B0329	B0429	B0529	B0629	B0729	B0829	B0929	B1029	 CN301

IBP Cables

Monitor	P/N	Conn.
<b>Drager-Siemens</b>		
5731281	DMS3042-04	
<b>Spacelabs</b>		
700-0028-00	DMS3012-16	
<b>Mindray</b>		
040-001029-00	DMS3013-12	
<b>GE</b>		
2005772-001	DMS3072-07	

IBP

## Temperature Probes

- Compatible with most of monitors
- Easy to use
- Latex free

### T 2 3 07 :

T: Temperature probes

1: Adult Skin type 2:Adult General/Rectal type 3: Pediatric skin type  
4:Pediatric general/rectal type)

3: 3.0m

07: YSI (01: Air shield; 02: Artema)

Warranty: one-year warranty for reusable temperature probes      Package: 1pcs/bag



## Temperature Probes

Monitor	Adult Skin	Adult General purpose	Pediatric Skin	Pediatric General purpose	Len.	Connector
<b>Yonker (E8)</b>	T1366	T2366	T3366	T4366	3.0m	 CN440
<b>Infinium DRE, OMNI</b>	T1367	T2367	T3367	T4367	3.0m	 CN314C
<b>Waltvick</b>	T1370	T2370	T3370	T4370	3.0m	 CN440
<b>Bionet</b>	T1371	T2371	T3371	T4371	3.0m	 CP2
<b>Edan</b>	T1372	T2372	T3372	T4372	3.0m	 CN406
<b>Air Shield</b>	T1301	T2301	T3301	T4301	3.0m	 CN401
<b>Artema-S&amp;W</b>	T1302	T2302	T3302	T4302	3.0m	 CN402
<b>GE-Marquette</b>	T1304	T2304	T3304	T4304	3.0m	 CN404
<b>Philips-HP</b>	T1305	T2305	T3305	T4305	3.0m	 CN405
<b>Creative</b>	T1310	T2310	T3310	T4310	3.0m	 CP2

TEMP

Temperature Probes

Monitor	Adult Skin	Adult General purpose	Pediatric Skin	Pediatric General purpose	Len.	Connector
Drager-Siemens	T1303	T2303	T3303	T4303	3.0m	 CN403
	T1322-D		T3322-D			 CN421
	T1322	T2322	T3322	T4322		 CN262
Drager C2000 VERSALET	T1303C	T2303C	T3303C	T4303C	3.0m	 CN456
Mindray	T1306	T2306	T3306	T4306	3.0m	 CN406
	T1318	T2318	T3318	T4318		 CN411
YSI 400	T1307	T2307	T3307	T4307	3.0m	 CN407
	T1307C	T2307C	T3307C	T4307C	3.0m	 CN407B
YSI 700	T1308	T2308	T3308	T4308	3.0m	 CN401
10k YSI	T1309	T2309	T3309	T4309	3.0m	 CN407

Temperature Probes


Monitor	Adult Skin	Adult General purpose	Pediatric Skin	Pediatric General purpose	Len.	Connector
Spacelabs	T1311	T2311	T3311	T4311	3.0m	 CN409
Ohmeda	T1312	T2312	T3312	T4312	3.0m	 CN410
David incubator HKN-90	T1317	T2317	T3317	T4317	3.0m	 CP6
David warmer (Ti2000)	T1360	T2360	T3360	T4360	3.0m	 CP4
	T1360B	T2360B	T3360B	T4360B		 CP2
Biolight	T1324	T2324	T3324	T4324	3.0m	 CN414
ATOM	T1315	T2315	T3315	T4315	3.0m	 CN05B
MEDIX	T1351	T2351	T3351	T4351	3.0m	 CN401
SAPS	T1316	T2316	T3316	T4316	3.0m	 CN407
Ohmeda 10K	T1320	T2320	T3320	T4320	3.0m	 CN424

TEMP








Temperature Probes

Monitor	Adult Skin	Adult General purpose	Pediatric Skin	Pediatric General purpose	Len.	Connector
DC 3.5mm plug	T1352	T2352	T3352	T4352	3.0m	 CN424
Medica D	T1353	T2353	T3353	T4353	3.0m	 CN401
Tary Mfg	T1354	T2354	T3354	T4354	3.0m	 CN412
Fanem	T1355	T2355	T3355	T4355	3.0m	 CN413
	T1155	T2155	T3155	T4155	1.5m	
	T1350	T2350	T3350	T4350	3.0m	 CN401
Disonmed	T1356	T2356	T3356	T4356	3.0m	 CP2
Comen	T1362	T2362	T3362	T4362	3.0m	 CN416
Zoncare	T1368	T2368	T3368	T4368	3.0m	 CN462
Contec	T1375	T2375	T3375	T4375	3.0m	 CN454

Temperature Probes

Monitor	Adult Skin	Adult General purpose	Pediatric Skin	Pediatric General purpose	Len.	Connector
Bistos	T1373	T2373	T3373	T4373	3.0m	 CN449

Fisher & Paykel Flow Sensor & Accessories




Monitor Brand	Compatible	P/N	Description	Photo	Original P/N
Fisher & Paykel	MR869 MR850	DMS9001	Temperature/Flow probe adaptor for use with circuits, 6pin 40degree		
	Mr850	DMS9002	Heater wire adaptor for dual-heated breathing circuits 900MR805		
	MR850	DMS9003	Heater wire adaptor for inspiratory heated breathing circuits 900MR806		
	MR700/730	DMS9004	Humidifier heater adapter cable, 4pin 2slot 80 degree, for disposable breathing circuit		
	MR730	DMS9005	Temperature flow sensor		
	900iw001	DMS9006	Pediatric skin temperature probe		
	MR850	DMS9007	Humidifier Single heater wire adaptor, 6pin 80degree		

TEMP

### Disposable Temperature Probes

Monitor Brand	Compatible Model	MEDKE P/N	Description	Monitor side	Probe side
Philips	All Philips temperature monitors that accept two-prong plugs	DMS5105 21091A	Disposable Skin-surface Probe,0.8m, 25°C@2.25KΩ	 CN424	
		DMS6105	Adult Esophageal/Rectal Probe,0.8m, 25°C@2.25KΩ		
YSI 400/700 Series	Mindray/Philips/Comen All China OEM brands	DMS5106 4499	Disposable Skin-surface Probe,0.8m, 25°C@2.25KΩ	 CN463	
		DMS6106	Adult Esophageal/Rectal Probe, 0.8m, 25°C@2.25KΩ		
Drager	Infant Incubator/ Radiant Warmer	DMS5122	Disposable Skin-surface Probe,0.8m, 25°C@2.25KΩ	 CN262	
		DMS6122	Adult Esophageal/Rectal Probe, 0.8m, 25°C@2.25KΩ		
	Infant Incubator/ Radiant Warmer	DMS5122-C	Disposable Skin-surface Probe,0.8m, 25°C@2.25KΩ	 CN262C	
		DMS6122-C	Adult Esophageal/Rectal Probe, 0.8m, 25°C@2.25KΩ		
Ohmeda	Giraffe, Omnibed, Panda iRes	DMS5112	Disposable Skin-surface Probe,0.8m, 25°C@10KΩ	 CN410	
		DMS6112	Adult Esophageal/Rectal Probe, 0.8m, 25°C@10KΩ		














### Smart Temp. Probe

Monitor Brand	Compatible	MEDKE P/N	Description	Photo	Original P/N
GE	GE Dinamap series V100	DMS7125	TurboTemp Oral Probe Blue		2008774-001
		DMS7125B	TurboTemp Oral Probe , Rectal ,Red		2008775-001
	Pro Services Pro1000	DMS7163	TurboTemp Oral Probe Blue		88012
		DMS7163B	TurboTemp Oral Probe Rectal ,Red		88016
WelchAllyn	WelchAllyn	DMS7164	TurboTemp Oral Probe Blue		








TEMP



Temperature Adapter Cable

Monitor Manufacturers	MEDKE P/N	Monitor side conn	Probe side conn	Len.
Philips-HP	DMS0205	 CN405	 CN423	2.4m
	DMS0005		 CN422	0.3m
GE-Marquette	DMS0204		 CN423	2.4m
	DMS0004-1	 CN404	 CN422	0.3m
	DMS0004-2		 CN422	0.3m
Mindray	DMS0206	 CN406	 CN423	2.4m
	DMS0006-1		 CN422	0.3m
Drager-Siemens	DMS0203	 CN403	 CN423	2.4m
	DMS0003		 CN422	0.3m

Temperature Adapter Cable

Monitor Manufacturers	MEDKE P/N	Monitor side conn	Probe side conn	Len.
Spacelabs	DMS0011-1	 CN409	 CN422	0.3m
	DMS0011-2		 CN422	0.3m
YSI400	DMS0207	 CN407	 CN423	2.4m
Artema-S&W	DMS0002	 CN402	 CN422	0.3m

TEMP



## Fetal Probes

- Compatible perfectly with famous fetal monitor
- Easy to insert and clear
- Latex free

**Warranty:** one-year warranty  
**Package:** 1pcs/bag



Monitor manufacturers	Connector	Options	P/N	OEM #	Compatible
OXFORD SONICAID TEAM		US Transducer (1.5Mhz)	DMS-001	8400-6919	TEAM/FM800
		US Transducer (2Mhz)	DMS-054	8400-6920	TEAM/FM800
		TOCO Transducer	DMS-016	8400-6921	TEAM/FM800
		Event marker	DMS-087		
Philips-HP		US Transducer	DMS-002	M1356A	1350/50A Series
		Event marker	DMS-003		50A series
		TOCO Transducer	DMS-012	M1355A	1350/50A Series
		Event marker	DMS-024	989803143411	Avalon FM series
		Philips Transducer Repair Cable	DMS0014		M2736A
Edan		US Transducer	DMS-011		CADENCE I
		US Transducer	DMS-007	MS3-109301(B)	CADENCE II F2/F6/F9/F3 SONATINA
		US Transducer	DMS-042		F2 (new)

Monitor manufacturers	Connector	Options	P/N	OEM #	Compatible
<b>Edan</b>		TOCO Transducer	DMS-006	MS3-31527-A1	F2/F3/F6/F9
		TOCO Transducer	DMS-005	TP01-RQ-1D	F9 express
		TOCO Transducer	DMS-005B		F9 (2017)
		Event marker	DMS-084		
<b>Bionet</b>		US Transducer	DMS-026		FC700
		TOCO Transducer	DMS-027		FC700
<b>Bistos</b>		US Transducer	DMS-031		BT300/MIDI
		TOCO Transducer	DMS0032		BT300/MIDI
<b>Sunray</b>		US Transducer	DMS-028		SRF618B5
		TOCO Transducer	DMS-029		SRF618B5
<b>Biocare</b>		US Transducer	DMS-020		FM-801
		TOCO Transducer	DMS-038		FM-801

Monitor manufacturers	Connector	Options	P/N	OEM #	Compatible
<b>Goldway</b>		US Transducer	DMS-025		CTG7 US probe
		US Transducer	DMS-033		UT3000A
		TOCO Transducer	DMS-035		G7 UT3000A/CTG7
<b>GE Corometric</b>		US Transducer	DMS-013	5700HAX	All series
		TOCO Transducer	DMS-039	2264HAX	All series
<b>Comen</b>		US Transducer	DMS-046		C20
		TOCO Transducer	DMS-048		STAR5000F
<b>Vconmin</b>		US Transducer	DMS-044		FM-321
		TOCO Transducer	DMS-045		FM-321
<b>Aeon</b>		US Transducer	DMS-040		A100P
		TOCO Transducer	DMS-041		A100P

Remark: please advise us brand ,model or connector image in order to match you correct transducer .as some OEM manufacturer have different version for same model



ESU Accessories

ESU Accessories

Photo	Description	P/N
	Reusable Bipolar Cable highly flexible silicone coated cable Autoclavable Euro Type to Banana 4.0mm	CP1002
	Reusable Bipolar Forceps cables ,3M highly flexible silicone coated cable Autoclavable USA 2-Pin Plug to 2*Banana 4.0mm	DMS1003
	Reusable silicone Monopolar Endoscopic Cable,3M highly flexible silicone coated cable Autoclavable 4mm Banana Plug > Female Banana	DMS1009
	Laparoscopy instrument adapter cable 3m,Autoclavable	DMS1015
	Aesculap reusable bipolar cableHIFI	DMS1009
	Reusable Bipolar Cable,3m Autoclavable EU fitting	DMS1015
	Wolf reusable Bipolar cable,3m Female Banana 4.0mm	DMS1019
	Grounding Pad Cable, 3M ,Hi-Fi plug or Erbe plug 6.3 Audio Plug with a single pin	DMS1004

Photo	Description	P/N
	Grounding Pad Cable , double 4.0 banana plug with a single pin,3M	DMS1005
	Grounding Pad Cable, 3M, REM plug "square" shape connector (2*φ2.35 female connectors with a single pin)	DMS1006
	Grounding Pad Cable, 3M "8" shape connector (2*φ2.35 female connectors with a single pin)	DMS1006B
	Disposable adult Monopolar Grounding pad	DMS1007
	Disposable adult Bipolar Grounding pad	DMS1008
	Reusable silicone patient plate with cable (HIFI), monopolar	DMS1012B
	Reusable silicone patient plate with cable (REM), monopolar	DMS1013B
	Disposable patient plate with cable (REM)	DMS1020

Photo	Description	P/N
	Single use Electrosurgical Pencil with the functions of cut and coagulation Button switch E.O sterilization	DMS1001
	Reusable use Electrosurgical Pencil with the functions of cut and coagulation Button switch E.O sterilization	DMS1001B
	Reusable Bipolar forceps, Straight 23CM, EU fitting	DMS10010
	Reusable Bipolar forceps, Straight 18CM, EU fitting	DMS1010B
	Reusable Bipolar forceps, Bayonet 24CM, EU fitting	DMS1011
	Reusable Bipolar forceps, Bayonet 20CM, EU fitting	DMS1011B
	Reusable Bipolar forceps, Bayonet 20CM, US fitting	DMS1014
	Reusable Bipolar forceps, Bayonet 24CM, US fitting	DMS1014B





Photo	Description	P/N
	<p>Nellcor® Veterinary Spo2 probe Sensor use for pet Ear clip or Tongue clip ,3M other compatible monitors are available such as Contec, BCI, Biolight .....</p>	<p>DMS3319A</p>
	<p>Veterinary ECG cable with flat Clip connection compatible with Contec C80Vet 3Lead/5Lead available AHA /IEC available</p>	<p>DMS3140V</p>
	<p>Veterinary two-in-one ECG and YSI400 temperature cable compatible with variety brands like Datex , Mindray,Philips , Comen</p>	<p>DMS3140A</p>

## Veterinary Accessories









Respiratory Consumables

Photo	Description	Band	P/N
<b>ETco2 Sensor</b>			
	External Sidestream CO2 Sensor	Respironics capnostat 5	DMS0039-B
		Philips	DMS0084-B
		Zoll	DMS0085-B
		Phasein (ISA AX+)	DMS0086-B
	Mainstream ETco2 Sensor Standard 8-pin 2-key LEMo connector (Bionet Bm5,Comen,Biolight)	Respironics capnostat 5	DMS0039
		Philips	DMS0084
		Zoll	DMS0085
		Phasein (IRMA CO2)	DMS0086
<b>Oxygen Sensor</b>			
	Original O2 sensor Oxygen cell for Draeger Savina 300	Dräger (Dräger Mx01049 )	111-000-0008
	Original O2 sensor Oxygen cell is a drop-in replacement for the 6803290, Narkomed G5, Cato, Narkomed 6000, Narkomed 6400, Anemone, Babylog 8000, Evita 2, Evita 2 Dura, Evita 4, Evita Series, Fabius, Fabius GS, Incubator 8000, Narkomed 1A, Narkomed 2, Narkomed 2B, Narkomed 3, Narkomed 4, Narkomed GS, Oxycom, Oxydig, Oxycom 100 D, PM8030, SA2 from Draeger	Dräger (Dräger 6850645 )	111-000-0009
	Compatible Oxygen Sensor 1.16 inch dia, Gold plated slip ring, 0 to 100%, White,Clear, 14 to 20 MV signal output, 700 to 1250 MBAR pressure, 12 sec response	Analytical Industries (PSR-11-915-2, Class R-23V, Unimed GO-02)	111-000-0014
	Compatible Oxygen Sensor Response 90%: 13 sec Expected Life: 60 months Traditional Style	Analytical Industries (PSR-11-917-M2, Class R-22Med)	111-000-0019

## Respiratory Consumables

Photo	Description	Band	P/N
<b>Oxygen Sensor</b>			
	Original O2 sensor HAMILTON-S1/G5 Galileo,Raphael	Hamilton	396200
	Original O2 sensor For gas supply in anesthesia and respiratory therapy (FiO2) incubators lung function analysis and in sports medicine	Envitec	DMS102-1
	Original O2 sensor Lead-free oxygen sensors for medical use Higher accuracy and reliability in response	Envitec OOM201	DMS001
	Original O2 sensor Lead-free oxygen sensors for medical use Higher accuracy and reliability in response	Envitec OOM202	DMS002
	Original O2 sensor Lead-free oxygen sensors for medical use Higher accuracy and reliability in response	Envitec OOM204	DMS003
	Original O2 sensor Lead-free oxygen sensors for medical use Higher accuracy and reliability in response	Drager (Drager 8606055 )	DMS004

## Respiratory Consumables

Photo	Description	Band	P/N
<b>Oxygen Sensor</b>			
	Original O2 sensor Maxtec Max-250b bare-wire oxygen sensor	Maxtec	DMS-250B
	Original O2 sensor Maxtec Max-250E replacement oxygen sensor for Maxtec MaxO2+AE, OM-25AE, OM-25ME, MaxBlend (Old Version), MaxVenturi, Carefusion Sentry Blender (Old Version)	Maxtec	DMS-250E
	Original O2 sensor Citytech CITY MOX1 Oxygen (O2) Gas Sensor	CITY MOX-1	DMS0X1
	Original O2 sensor Citytech CITY MOX2 Oxygen (O2) Gas Sensor	CITY MOX-2	DMS0X2
	Original O2 sensor Citytech CITY MOX3 Oxygen (O2) Gas Sensor	CITY MOX-3	DMS0X3
	Original O2 sensor Citytech CITY MOX4 Oxygen (O2) Gas Sensor	CITY MOX-4	DMS0X4



## Respiratory Consumables

Photo	Description	Brand	P/N
<b>Oxygen Sensor</b>			
	Original O2 sensor 3.5mm Mono Jack Connector, White/Clear	R-17MED	R-17MED
<b>Water Trap</b>			
	Waterlock 2 Trap Water (6872130)	Drager	DMS0047
	Dryline water trap, Adult/Pediatric	Mindray	DMS0075
	Dryline™ II water trap used with the 1x Sidestream CO2 module, Adult/Pediatric	Mindray	DMS0073
<b>Sampling Cannulas</b>			
	GE Oral/ Nasal CO2 Sampling O2 Delivery Cannula, Adult,3m	GE	73319-HEL
	CO2 nasal sample cannula	Mindray	DMS0038

## Respiratory Consumables

Photo	Description	Brand	P/N
<b>Breathing Circuit</b>			
	Disposable Breathing Circuit, Pediatric 1.5m single tube , Heating Type		DMS0064
	Adult Breathing Circuit, Silicone		DMS0054
	Pediatric Breathing Circuit, Silicone		DMS0054B
	Newborn Breathing Circuit, Silicone		DMS0054C
	Disposable Anesthesia Ventilator Breathing Circuit Extension Tube, Adult		DMS0068
<b>Ventilator Accessories</b>			
	flow sensor	GE	2089610-001-S

## Respiratory Consumables

Photo	Description	Brand	P/N
<b>Ventilator Accessories</b>			
	Original Drager 8403735 Spirolog Flow Sensor 5pcs/pouch	Drager	8403735
	Compatible Hamilton 281637 Adult Flow sensor	Hamilton	DMS0048
	Compatible Hamilton 260177 Infant Flow sensor	Hamilton	DMS0048B
	Fabian HFO Acutronic Flow Sensor Cable	Fabian	DMS0071
	Drager BabylogVN500 flow Sensor Cable	Drager	DMS0045 8416600
	neonatal Drager Babylog 8000 flow Sensor Cable	Drager	DMS0082 8409626

## Customer Service Guide

### 1. How to order

Call:

Dial 0086-755-23463462 during our business hour Monday to Friday 9:00 AM - 18:00 PM(Beijing Time).

Fax:

Fax your order to us: 0086-755-29553084 , Faxed orders are accepted 24hours a day.

Email:

Email your order to our sales or info@medke.com

Online

<http://stores.ebay.com/medke>

Notice:

Please inform us monitors model & name, Medke P/N and quantities while placing order

If you have any questions about the product, please contact with our sales departments

### 2 . Payment Methods

Wire Transfer

If you would like to remit payment directly to our bank, please transfer to the account exactly as Proforma invoice.

To expedite the payment process, please reference the Invoice number on the transfer documentation. And send your remitted paper or Bank slip to our sales via email when finished your payment. the following account is valid:

Bank Account (US\$ only) Information

Account No.: 4000 0551 0910 0071 650

Beneficiary Name: Shenzhen Medke Technology Co., Ltd

Bank Name: Industrial And Commercial BANK OF China, Shenzhen

Swift Code: ICBKCNBJSZN

Beneficiary Add.: RM.A, 4/F,Bldg.A1,Anle Ind. Zone, Hangcheng RD., Baoan Dist.,Shenzhen,China

Note: Customer is responsible for both bank charge and the full amount shown on invoice must be remitted as our offer is ExW shenzhen.

For small value of payment, please send it via western union or PayPal

For more choices, please contact with our sales service

### 3. Warranty Terms

No warranty for disposable products, unless any defective found during you first use

6months after delivery date for IBP cable, Reusable NIBP Cuff, EEG cable, Internal ECG cable, Suction/Limb clamp electrodes

one-year after delivery time for reusable Spo2 sensor &Ext. Cable, ECG/EKG cables and Leadwire, Temperature probes,

Fetal probes, Electrosurgical connecting cables

Original products refer to original manufacturer

### 4. Replace or Repair

All the returned goods should have S.N and P/N no matter it is for repair or replace

Freight Charges: Customer is responsible for all transport cost related to the return of items under warranty or for repair.

we cannot reimburse freight costs or other related transport costs for any returned goods. Freight collected shipments are not accepted.

★ If Medke determines that returned items are covered under warranty, then the items will be shipped back to our factory freight prepaid via DHL, TNT, UPS and Federal Express;

★ If items are not covered under warranty, customer is responsible for repairing costs and return shipping costs;

★ Merchandise returned due to customer order error will be charged for 10% (within 15 days) or 20% (within 90days) administration fee;

★ Medke will not accept returns after 6 months from the delivery date;

Notice: The listed products might be modified or discontinued without advising in advance, if you want to confirm, please contact with our sales

### 5. Statements

The brands and trademarks involved in the catalogue are to illustrate the compatibility of Medke products, and the brands or trademarks belong to the corresponding manufacturers. In order to avoid unnecessary disputes, as a customer of Medke, kindly please pay attention to whether the products sold in this area have product patent problems.